

## Mobile Fever Screening Station

### Powered Mobile Workstation with 24" All-in-One PC and Thermographic Turret Body Temperature Measurement Camera.

An all-in-one powered mobile workstation providing rapid, precise and contactless measurement of human body temperature using facial recognition and thermal imaging technology to provide fever screening for public locations with multiple person entry. This mobile workstation provides a complete standalone and 'ready to roll' solution, without any need for installation. Multiple adjustment features enable highly flexible positioning options.

#### Applications:

- Hospitals
- Airports
- Train stations
- Sports venues
- Theatres
- Shopping centres
- Offices

#### Key features:

##### Cart:

- Height-adjustable Workstation on Wheels (WOW) with multi-adjustable keyboard and mouse tray
- Lithium-ION battery and power inverter for 6-8 hrs run time
- LCD battery indicator display

##### Camera:

- 160 × 120 resolution in thermal channel, high sensitivity sensor
- Temperature range: 30°C to 45°C; temperature accuracy:  $\pm 0.5$  °C
- Reliable temperature exception alarm function
- 3D DNR, image detail enhancement
- Audio alarm.

Patent pending.



## Mobile Fever Screening Station

Powered Mobile Workstation with 24" All-in-One PC and Thermographic Turret Body Temperature Measurement Camera.



Height adjustable and ergonomic



Storage basket



Hand Sanitiser holder

Patent pending.