

Better **solutions.**
Better **care.**

Agile Medical

Mobile Theatre Equipment



Welcome

Agile Medical is dedicated to improving workflow and efficiency in healthcare environments through the implementation of mobile solutions that have been developed with Ergonomics, Health & Safety and Infection Control compliance in mind.

Contents

Narrow Anaesthetic Trolley	4	Surgical Steps	19
LX Anaesthetic Trolley	5	Bowl Stands	20
Double Anaesthetic Trolley	6	Sterile Storage Carts	21
Theatre Recovery Trolley	7	Infusion Stands	22
Procedure Trolley	8	Mobile Theatre Cabinet	23
Difficult Airway Trolley	9	Standard Linen Trolleys	24
Monitor Arms	10	Mobile Theatre Trolley	25
Heavy Duty Roll Stand	12	Theatre Stools	26
White Instrument Trolley	13	Theatre Operator Stools	27
Stainless Steel Instrument Trolleys	14	Mobile Privacy Screens	28
Special Instrument Trolley	16	Theatre Trolley Accessories	29
Special Theatre Trolley	17	Theatre Drawer Accessories	32
Mayo Tables	18	NHS How to Order	33

Agile Medical
Better **solutions.**
Better **care.**

Narrow Anaesthetic Trolley

Mobilising the storage of anaesthetic equipment reduces staff downtime and prevents ‘dead stock’ building up. The narrow anaesthetic trolley provides these benefits in a compact, space optimising solution that can be configured to your specific requirements.

- Security and privacy
- Staff safety in mind
- Infection prevention - easy clean
- Anti bacterial powder / nylon coating
- 5 year guarantee



Supplied accessories:

Label Holders



Push handle



2 x Divider sets



Slide pull-out extension table



Trolley Options



Code **NC3405_ (1 or 2)**
Dimensions: 520w, 540d, 946h (mm)
5 drawer layout
(2x 3", 2x 6", 1x 9")

- Colours:
- 1 White
 - 2 Light Blue

For all accessories see pages 29-32.



Code **NC34GAR_ (1 or 2)**
Dimensions: 520w, 540d, 946h (mm)
5 drawer layout
(2x 3", 2x 6", 1x 9")

- Colours:
- White Paediatric (Garden Theme)
 - Light Blue Paediatric (Garden Theme)

For all accessories see pages 29-32.

LX Anaesthetic Trolley

Tried and tested in the theatre departments of many leading NHS hospitals, this range exceeds the infection control and durability requirements of a theatre environment.

- Security and privacy
- Staff safety in mind
- Infection prevention - easy clean
- Anti bacterial powder / nylon coating
- 5 year guarantee



Supplied accessories:

Label Holders



Push handle



2 x Divider sets



Slide pull-out extension table



Trolley Options



Code **LX31ANE_ (1 or 2)**
Dimensions: 690w, 540d, 930h (mm)
4 drawer layout
(1x 3", 2x 6", 1x 9")

- Colours:
- 1 White
 - 2 Light Blue

For all accessories see pages 29-32.



Code **LX34ANE_ (1, 2 or 3)**
Dimensions: 690w, 540d, 1010h (mm)
5 drawer layout
(2x 3", 2x 6", 1x 9")

- Colours:
- 1 White
 - 2 Light Blue
 - 3 Blue

For all accessories see pages 29-32.



Code **LX34GAR_ (1 or 2)**
Dimensions: 690w, 540d, 1010h (mm)
5 drawer layout
(2x 3", 2x 6", 1x 9")

- Colours:
- White Paediatric (Garden Theme)
 - Light Blue Paediatric (Garden Theme)

For all accessories see pages 29-32.

Double Anaesthetic Trolley

The optimum storage solution for all your anaesthetic equipment. The combination of LX and NC storage on one trolley brings greater flexibility with equipment storage and ensures all the vital equipment is right at hand throughout the entirety of the procedure.

- Security and privacy
- Staff safety in mind
- Infection prevention - easy clean
- Anti bacterial powder / nylon coating
- 5 year guarantee



Supplied accessories:

Slide pull-out extension table

10 x Label Holders

4 x Divider sets

Push handle

Trolley Options



Code **DC34092**
Dimensions: 1125w, 545d, 1010h (mm)
Colours: **2** Light Blue
For all accessories see pages 29-32.



Code **DC34GAR_ (1 or 2)**
Dimensions: 1125w, 545d, 1010h (mm)
Colours: White Paediatric (Garden Theme)
 Light Blue Paediatric (Garden Theme)
For all accessories see pages 29-32.

Theatre Recovery Trolley

The Theatre Recovery Trolley is designed to assist healthcare professional staff in Recovery Departments – Supplied c/w 100% telescopic pull-out keyboard shelf, screen holder and server/box holder.

- Security and privacy
- Staff safety in mind
- Infection prevention - easy clean
- Anti bacterial powder / nylon coating
- 5 year guarantee



Supplied accessories:

4 x Label Holders

2 x Divider sets

PC holder

4 power socket and lead

Slide pull-out extension table

Support holder & computer arm

Trolley Options



Code **LX34TRT_ (2 or 3)**
Dimensions: 690w, 540d, 1010h (mm)
5 drawer layout (2x 3", 2x 6", 1x 9")
Colours: **2** Light Blue
3 Blue
For all accessories see pages 29-32.

Procedure Trolley

Our Procedure Trolleys can be arranged to suit a multitude of medical departments and situations.

- Security and privacy
- Staff safety in mind
- Infection prevention - easy clean
- Anti bacterial powder / nylon coating
- 5 year guarantee



Supplied accessories:

Trolley Options

Code

LX28PRO_ (1 or 2)

Dimensions:

690w, 540d, 850h (mm)
4 drawer layout
(1x 3", 3x 6")

Colours:

1

White

2

Light Blue

For all accessories see pages 29-32.

Code

LX31PRO_ (1 or 2)

Dimensions:

690w, 540d, 930h (mm)
4 drawer layout
(1x 3", 2x 6", 1x 9")

Colours:

1

White

2

Light Blue

For all accessories see pages 29-32.

Code

LX34PRO_ (1 or 2)

Dimensions:

690w, 540d, 1010h (mm)
5 drawer layout
(2x 3", 2x 6", 1x 9")

Colours:

1

White

2

Light Blue

For all accessories see pages 29-32.

Difficult Airway Trolley

The Difficult Airway Trolley can be customised, with 4 or 5 drawer option available and a zero risk emergency breakaway locking system and many other unique features. Each trolley is supplied with drawer front labels. These labels help clinical and non-clinical staff identify the correct procedure in a difficult airway or intubation situation.

- Staff safety in mind
- Infection prevention - easy clean
- Anti bacterial powder / nylon coating
- 5 year guarantee



Supplied accessories:

Trolley Options

Code

LX28AIR_ (5 or 8)

Dimensions:

690w, 540d, 850h (mm)
4 drawer layout
(1x 3", 3x 6")

Colours:

5

Red

8

Dark Green

For all accessories see pages 29-32.

Code

LX31AIR_ (5 or 8)

Dimensions:

690w, 540d, 930h (mm)
4 drawer layout
(4x 6")

Colours:

5

Red

8

Dark Green

For all accessories see pages 29-32.

Code

LX34AIR_ (5 or 8)

Dimensions:

690w, 540d, 1010h (mm)
4 drawer layout
(3x 6", 1x 9")

Colours:

5

Red

8


Dark Green

For all accessories see pages 29-32.

Monitor Arms

Mounting solutions that set new standards.

Compliant with the highest medical regulations, plus integrated cable routing and a 5-year maintenance-free warranty. What's really at our core is Infection Prevention. CIM arms are clinically proven to be the best in the world for Infection Prevention, according to independent test results.



Infection prevention through easy-to-clean surfaces.

Clinical test results, comparing all of the leading brands in medical mounting solutions, revealed CIM arms to be at least four times better for advanced hygiene control.



Integrated cable management

We meet the growing demands of hygiene and aesthetics in the medical environment while protecting cables from damage and accidental disconnection.



Silky ergonomics

Smooth adjustment and maneuverability with perfect placement of the monitors and devices.



Flexible positioning

Large swiveling and turning ranges for optimum positioning of the devices.



Space saving

Optimal parking positions through flexibility.



Zero maintenance

Guaranteed maintenance-free operation during the standard 5-year warranty.



6-fold safety factor

The highest safety standard on the market for patient, user and device safety.



TÜV-certified safety

Independently tested systems that comply with current EN 60601 regulations.



Precise engineering

Guarantees a level screen in all viewing positions. Finesse at your fingertips.



CIM mounting solutions are compatible with patient monitors from a wide range of medical device manufacturers.

To discuss the various versions, colour schemes, and mounting options available, please get in touch with our product support team at sales@agilemedical.co.uk.

Heavy Duty Roll Stand

The Heavy-Duty Roll Stand can be configured according to your needs with an extensive range of accessories available to match your requirements. All individual solutions are tested and certified before delivery and are tip tested in conformance with EN 60601-1, 3rd Edition.



Accepts:

- Drawers and shelves:** Single and double drawer consoles or simply shelves to place your equipment.
- Support arms:** CIM height adjustable and non-height adjustable arms to mount monitors, devices and other medical equipment.
- Accessories:** A huge variety of accessories is available: cable hooks, baskets, IV poles.

Monitor Presentations

Which monitor would you like to mount?
Agile Medical offers solutions for a wide range of medical device manufacturers.



White Instrument Trolley

Good quality Instrument Trolleys for hospitals and health care centres. These Trolleys are available in Single or Double width.

- Staff safety in mind • Infection prevention - easy clean • Anti bacterial powder / nylon coating • 5 year guarantee



Trolley Options



Code

IT3502SW

Dimensions:

500w, 500d, 900h (mm)



Code

IT3504DW

Dimensions:

1000w, 500d, 900h (mm)



Single ABS Top

Code

IT5050T

Dimensions:

500w, 500d (mm)



Double ABS Top

Code

IT1050T

Dimensions:

1000w, 500d (mm)

Stainless Steel Instrument Trolley

Instrument Trolley – Strong (grade 304) stainless steel seamless welded construction with rounded corners (no sharp edges), swivel anti-static castors and rubber buffers.

- Staff safety in mind
- Infection prevention - easy clean
- 5 year guarantee



Features:

Anti-static castors

Flush-welded, smooth finish and easy clean

Buffers

Trolley Options

Code **SS4545F**
Dimensions: 450w, 450d, 840h (mm)

Code **SS7545F**
Dimensions: 750w, 450d, 840h (mm)

Code **SS9145F**
Dimensions: 910w, 450d, 840h (mm)

Stainless Steel Instrument Trolley

Instrument Trolley – Strong (grade 304) stainless steel seamless welded construction. Flush welded shelves with raised rim containment edges, swivel anti-static castors and rubber buffers.

- Staff safety in mind
- Infection prevention - easy clean
- 5 year guarantee



Features:

Anti-static castors

Flush welded shelves

Buffers

Special Instrument Trolley

Tall Instrument Trolley – Strong (grade 304) stainless steel seamless welded construction with rounded corners (no sharp edges), swivel / lockable castors. Available for 900mm and 1000mm workshop heights.

- Staff safety in mind
- Infection prevention - easy clean
- 5 year guarantee



Trolley Options



Code **SIT0994**
Dimensions: 900w, 450d, 900h (mm)



Code **SIT1094**
Dimensions: 900w, 450d, 1000h (mm)

Special Theatre Trolley

Special Theatre Trolley – Strong (grade 304) stainless steel seamless welded construction. All shelves with raised rim containment edges, swivel / lockable castors and removable 15L stainless steel buckets (available separately).

- Staff safety in mind
- Infection prevention - easy clean
- 5 year guarantee



Options



Special Theatre Trolley
Code **CIT0947**
Dimensions: 770w, 450d, 900h (mm)



15L Bucket
Code **NSB300**
Diameter: 295mm
Height: 310mm

Mayo Tables

Easy to use, mobile Mayo Tables providing a place for sterile instruments and supplies used during surgery. Fully stainless-steel welded construction providing compatibility with all standard hospital cleaning products.

- Staff safety in mind
- Infection prevention - easy clean
- 5 year guarantee



Options



Mayo Table with manual height adjustment
Code **HMTD46**

- Stainless steel Frame (AISI 304)
- ST/ST top with perimetrical rim (50cm x 70cm WxD)
- Manual height adjustment (80-120 cm)
- 4 medical grade castors (2 braked)



Mayo Table with hydraulic height adjustment
Code **HMTD48**

- Stainless steel Frame (AISI 304)
- Removable ST/ST top with perimetrical rim (50cm x 70cm WxD)
- Hydraulic foot pedal height adjustment (80-120 cm)
- 4 medical grade castors (2 braked)



Mayo Table with hydraulic height adjustment
Code **HMTD50**

- Stainless steel Frame (AISI 304)
- ST/ST top with perimetrical rim (50cm x 70cm WxD)
- Hydraulic foot pedal height adjustment (80-120 cm)
- 4 medical grade castors (2 braked)


Surgical Steps

Surgeon Steps – Stable 4-footed single step construction, grade 304 stainless steel seamless welded frame construction with removable stainless steel ribbed tread for easy cleaning. 2 non-slip, non-marking adjustable rubber feet and 2 conductive polymer screw-in feet.

- Staff safety in mind
- Infection prevention - easy clean
- 5 year guarantee



Features:



Options



Code **STE1631**
Dimensions: 600w, 350d, 130h (mm)



Code **STE1632**
Dimensions: 600w, 350d, 220h (mm)



Code **STE1633**
Dimensions: 600w, 350d, 300h (mm)



Code **STE26624**
Dimensions: 600w, 650d, 220/440h (mm)

Bowl Stands

Manufactured from grade 304 stainless steel, our bowl stands and kickabouts exceed infection control requirements for within sterile environments. All parts can easily be cleaned and sanitised.

- Staff safety in mind
- Infection prevention - easy clean
- 5 year guarantee



Options



Single Bowl Stand
Code **AM/BSD-79**

- Dimensions: 350w, 350d, 800h (mm)
- Stainless steel Frame (AISI 304)
 - 5-star ST/ST base with swivel castors (2 braked)
 - 1 removable ST/ST bowl



Double Bowl Stand
Code **AM/BSD-80**

- Dimensions: 700w, 350d, 800h (mm)
- Stainless steel Frame (AISI 304)
 - 5-star ST/ST base with swivel castors (2 braked)
 - 2 removable ST/ST bowls



Single Kick Bucket
Code **AM/KBM-12**

- Dimensions: 750 Ø, 335h (mm)
- ST/ST Base and Bucket (AISI 304)
 - 5 antistatic castors, Ø 75 mm (2 braked)
 - 1 removable 12L bucket
 - Available with and without handle

Sterile Storage Carts

Single & Double Sided Storage Carts – Strong (grade 304) stainless steel seamless welded frame with hook configuration to accommodate 3 – 4 baskets, hooks tested to 50kg load capacity, swivel / lockable castors.

- Staff safety in mind
- Infection prevention - easy clean
- 5 year guarantee



Options



Single Storage Cart
Code **STC300**

- Dimensions: 610w, 470d, 1450h (mm)



Double Storage Cart
Code **STC600**

- Dimensions: 610w, 680d, 1450h (mm)

Infusion Stands

Light, medium and heavy duty Infusion Stands to suit a variety of applications
– Robust design, stable 5-star base, easy to clean and simple to operate.

- Staff safety in mind
- Infection prevention - easy clean
- 5 year guarantee



Options



Light Duty Infusion Stand
2 Hooks **NS17230**
4 Hooks **NS17229**
– Colour of Plastic Parts: Blue
– Material: Stainless Steel
– Height Adjustment: 1350-2150mm
– Base Diameter: 635mm
– Base Weight: 3.5kg



Medium Duty Infusion Stand
2 Hooks **ICS17230**
4 Hooks **ICS17229**
– Colour of Plastic Parts: Blue
– Material: Stainless Steel
– Height Adjustment: 1500-2350mm
– Base Diameter: 680mm
– Base Weight: 8.5kg / 6.5kg



Heavy Duty Infusion Stand
Code **PCI172954**
– Specifically for pump attachments
– 4 Hooks – Stainless Steel
– Colour of Plastic Parts: Black
– 5 star 630mm diameter base
– Height Adjustment: 1430-2260mm (soft drop system)
– Load capacity: 32kg
– Twin castors - electrically conductive and lockable

Mobile Theatre Cabinet

Single or Double HTM71 Mobile Base Cabinets are supplied with lockable doors and internals. Specially designed to fit below 900mm worktop heights.



Cabinet Options



Single B40 Cabinet
Code **SMC846_**
Dimensions: 450w, 650d, 860h (mm)
Comes with: - 1 x Door (lockable)
- 4 x 600x400x100mm Poly Trays
- 1 x 600x400x200mm Poly Basket



Double B40 Cabinet
Code **DMC896_**
Dimensions: 900w, 650d, 860h (mm)
Comes with: - 2 x Doors (lockable)
- 8 x 600x400x100mm Poly Trays
- 2 x 600x400x200mm Poly Basket

The above cabinets are suited for 600mm insertions.
Also available to suit 400mm insertions.

Type A - 600mm insertion



Type B - 400mm insertion



Various door colours available (contact us for full list):
To order add colour code required (e.g; DMC896-U522)

W980	U708	F501	U522	U539	M047	U400	U311
White	Light Grey	Titanium	Horizon Blue	Steel Blue	Mousse	Powder Violet	Burgundy

Standard Linen Trolleys

A range of easy-to-clean, mobile trolleys to provide enclosed storage of clean linen. One of our best selling ranges for Endoscopy, Theatres and Day Surgery.

- Staff safety in mind
- Infection prevention - easy clean
- 5 year guarantee



Trolley Options



TYPE A
Code **LT1085**
Dimensions: 740w, 535d, 1000h (mm)
Comes with 1 central shelf.



TYPE B
Code **LT1285**
Dimensions: 740w, 535d, 1190h (mm)
Comes with 2 shelves.



TYPE C
Code **LT1485**
Dimensions: 740w, 535d, 1400h (mm)
Comes with 2 shelves.

Double Linen Cabinet

Code **LTC875**
Dimensions: 700w, 540d, 800h (mm)
Comes with two internal compartments with adjustable shelves.



Mobile Theatre Trolley

An attractive and specially designed HTM71 mobile cabinet for the theatre departments. High impact strength, easy-to-clean and resistant to chemicals and high energy radiation (gamma and x-rays).



Various door colours available:
To order add colour code required (e.g; MTC1675-6)

1	2	3	4	5	6
Yellow	Orange	Red	Teal	Blue	Grey

Trolley Options



Mobile Theatre Trolley
Code **MTC1675_**
Dimensions: 740w, 580d, 1655h (mm)
Body: Stainless Steel
Door Frame: PE500 with part transparent
Capacity: - 9 x 600x400x100mm Trays
or
- 8 x 600x400x100mm Trays
- 1 x 600x400x200mm Tray

Theatre Stools

Medical Line range with comfortable 40cm round seat, attractive colour options and 10cm height adjustment creating an optimum working area of between 80-90cm.



Options

Medical Stool (Round Seat)
Code **MR6101_**

Height adj. 570-760mm
Seat diameter 400mm

Medical Stool (Round Seat) with backrest
Fixed seat angle **MR6211F_**
Adjustable seat angle **MR6211T_**

Height adj. 570-760mm
Seat diameter 400mm

Medical Stool (Round Seat) with backrest
Fixed seat angle **MR6221F_**
Adjustable seat angle **MR6221T_**

Height adj. 570-760mm
Seat diameter 400mm

Height Adjustment by Foot
Code **FA6000**



Various colour fabrics available.
To order add colour code required (e.g; MR6211F-**K89**)

K89

K05

K85

K195

K80

K184

K81

K93

K87

K170

Dark Blue

Black

Anthracite

Grey

Light Blue

Turquoise

Bright Green

Bright Orange

Bright Red

Burgundy Red

All models are also available with the **Ergo seat**.
Contact us for product details.



Theatre Operator Stools

Our range of operator stools designed and developed with input from leading ergonomists and surgeons to support theatre operators throughout procedures, by creating a stable and ergonomic sitting posture.



Options

Amazon Saddle Stool
Fixed seat angle **MAS6311F_**
Adjustable seat angle **MAS6311T_**

Height adj. 590-770mm
Seat appr. 340w, 300d (mm)

Jumper Saddle Stool
Fixed seat angle **MJS6311F_**
Adjustable seat angle **MJS6311T_**

Height adj. 590-770mm
Seat appr. 440w, 300d (mm)

Ergo Support Surgeons' Stool
Code **ME6331TFHB**

- Adjustable seat angle
- Foot operated height adj. 520-710mm
- Central foot brake
- 2 swivel, adjustable armrests

Height Adjustment by Foot
Code **FA6000**



Various colour fabrics available.
To order add colour code required (e.g; MAS6311-**K89**)

K89

K05

K85

K195

K80

K184

K81

K93

K87

Dark Blue

Black

Anthracite

Grey

Light Blue

Turquoise

Bright Green

Bright Orange

Bright Red

Mobile Privacy Screens

Quick and easy to open out and close away, these mobile privacy screens are ideal for having ready-to-hand in a multi-bay recovery area to help maintain the privacy of patients.



Options



3 Panel Privacy Screen
Code **AM31400**
Dimensions: 1650L x 1680H (mm)
– 3 white multi-use vinyl panels
– Aluminium tube frame
– Swivel Castors Ø50mm



4 Panel Privacy Screen
Code **AM31401**
Dimensions: 2220L x 1680H (mm)
– 4 white multi-use vinyl panels
– Aluminium tube frame
– Swivel Castors Ø50mm



Theatre Trolley Accessories



Anaesthesia Label Dispenser
Clear plastic anaesthesia label dispenser.
Dimensions: 300 x 127 x 140mm
Code: **AMDSP-HC10**



Horizontal 3 box glove holder – clear ABS
Suitable for narrow trolleys.
Dimensions: 260 x 395 x 105mm
Code: **GH2310W**



Vertical 3 box glove holder – clear ABS
Dimensions: 400 x 258 x 105mm
Code: **GH4210W**



IV pole
Adjustable IV pole
Code: **IVP100**



Tilt Bin Organiser (Half Size)
Suitable for narrow trolleys.
Dimensions: 300 x 111 x 230mm
Code: **TBS312**



Tilt bin organiser (small)
600 x 111 x 230mm
Code: **TBS612**



Tilt bin organiser (large)
600 x 161 x 168mm
Code: **TBL611**



Tilt bin organiser (extra large)
600 x 200 x 168mm
Code: **TBX612**



Accessory bridge
800mm
Code: **AB800**



Glove dispenser & back rail
160 x 260 x 95mm
Code: **GDB119**



Glove dispenser & side rail
160 x 260 x 95mm
Code: **GDS119**




3 box glove holder & side rail
400 x 260 x 85 xmm
Code: **GDS428**




Bulk storage bins (L) & back rail
75 x 105 x 190mm
Code: **BSB711**


Theatre Trolley Accessories




Waste bin & side rail
265 x 255 x 460mm
Code: **WBS224**




Waste container w/lid & side rail
255 x 170 x 315mm
Code: **WCS213**




Garbage bin & side rail
255 x 170 x 315mm
Code: **GBS213**




Oxygen holder (L) & side rail
125 x 125 x 700mm (E Size)
Code: **OTL117**




Suction Machine Holder & side rail
80 x 400 x 200mm
Code: **SMH300**




Catheter holder & side rail
100 x 100 x 480mm
Code: **CHS114**




Multi-storage & side rail
Code: **MSL100**




File holder & side rail
571 x 397 x 150mm
Code: **FHS531**




Sharps bin holder & side rail
Code: **SBH100**




Catheter holder B & side rail
Type B (Vertical Storage)
600 x 390 x 150mm (HWD)
Code: **301106/SR**




Apron roll dispenser & side rail
Can be trolley mounted
162 x 162 x 410mm
Code: **306171/SR**




Monitor arm and laptop holder
Code: **LAP100**




Sharps bin holder & side rail
100 x 190 x 110mm
Code: **SBH200**



Sharps bin holder
350 x 170mm
Code: **SBH300**




Gel holder
Code: **GMH100**




Monitor arm and screen holder
Code: **COM100**



Arm support holder
Support the laptop or computer arm and fix to the trolley
Code: **ASH100**




PC holder and 2 x side rails
Code: **PCH400**



Foldaway side drawer unit
ABS 3 side files (open or closed)
Code: **FAS416**



Stainless steel defibrillator shelf (extra strength)
330 - 410 x 330mm
Code: **DSS433**



Oxygen holder (S) & side rail
145 x 125 x 403mm (CD Size)
Code: **OTS114**




Side rails
450mm
Code: **ASR450**



Back rails
LX Trolleys - 580mm
Code: **ABR550**

DC Double Trolleys - 980mm
Code: **ABR980**



Trolley keys
Code: **LXK200**



Magnetic label holders
80 x 50 x T5mm
6pcs/pack
Code: **MLH855**

Theatre Drawer Accessories



Straight divider set
76mm fixed divider set in drawer (3 x short + 3 x long + 2 x strip)
Code: **LXF076**



Straight divider set
152mm fixed divider set in drawer (2 x short + 2 x long + 2 x strip)
Code: **LXF152**




Adjustable divider set
76mm adjustable divider set in drawer (3 x short + 6 x adjustable + 2 x 76mm strip + 6 x T clips)
Code: **LXA076**



Adjustable divider set
152mm adjustable divider set in drawer (3 x short + 4 x adjustable + 2 x 152mm strip + 4 x T clips)
Code: **LXA152**



Laptop/keyboard drawer (76mm)
All colours. Available for all LX trolleys (690mm wide)
Code: **PCB690**



ABS trays - small
Size: 71 x 570 x 396mm
Code: **LXT764**



ABS trays - large
Size: 152 x 570 x 396mm
Code: **LXT164**



Straight divider set (for trays)
76mm fixed divider set (3 x short + 3 x long)
Code: **FTD076**



Straight divider set (for trays)
152mm fixed divider set (2 x short + 2 x long)
Code: **FTD076**



Divider set for narrow trolleys
Size: 71 x 399 x 398 mm
Code: **NAA076**



Divider set for narrow trolleys
Size: 152 x 399 x 398 mm
Code: **NAA152**



How to Order:

If the NHS Trust you work at has Agile Medical set up as a supplier, you will just need a quotation from Agile Medical for the products you require.

Depending on how the Trust internal order process works, the quotation should be used to raise a requisition internally for Procurement to raise an official Purchase Order.



Orders can be sent to our **sales@agilemedical.co.uk** email address.



Or alternatively send as a hard copy in the post to:
Agile Medical, South Humberside Industrial Estate, Estate Rd 2, Grimsby, DN31 2TG



For orders over the phone, our representatives will require your full contact details, invoice and delivery address as well as the Purchase order number. A proforma invoice with our company bank details will be sent to your Finance team to make payment.



Payment can be made via Bank Transfer or cheque. Unfortunately at this time, we do not have the facility to accept payment over the phone.

Agile Medical

Head Office

Agile Medical
Estate Road 2
South Humberside Industrial Estate
Grimsby
DN31 2TG
UK

London Office

Unit 1 Kimberly Business Park
Blackness Lane
Keston
BR2 6HL
UK

Phone: +44 (0)1472 266 190

Email: sales@agilemedical.co.uk

Web: agilemedical.co.uk