

SterileFLAT®

Antibacterial Medical Keyboards & Mice

**ANTIBACTERIAL
BACKLIT
MEDICAL
KEYBOARD**



Reduce Cross Infections by upgrading to the SterileFLAT Antibacterial Backlit Medical Keyboard – SF20BL

SterileFLAT Antibacterial Backlit Medical Keyboards are an invaluable tool in reducing the risk of cross infection in any environment.

Shared keyboards on the wards or other clinical environments can become potential breeding grounds for Healthcare associated infections (HAI). Infections that are acquired in hospitals or as a result of healthcare interventions.

SterileFlat Keyboards offer a conventional style keyboard with Nano Silver impregnated antibacterial silicone latex free membrane barrier which can be used in all areas of clinical and medical environments.

They can be regularly washed with a recognised antibacterial or alcohol based cleaning agent every day, thus reducing the spread of MRSA, C. dif, Norovirus, E. coli, Salmonella and Airborne Viruses. Find out more at www.admor.co.uk

Key features/benefits of the SterileFLAT Keyboard:

- Conventional style keyboard with five levels of brightness.
- Very low profile keys with angled edges and wide caps to aid effective cleaning.
- Can be cleaned with Antibacterial and Alcohol based agents or Chlor-clean a chlorine based agent as used in many Hospitals.
- Fully sealed to IP-68 standard.
- Keyboard face Antibacterial. The membrane combines active Nano Silver technology with a flexible but strong latex free silicone barrier cover.
- On/Off power button.
- CE and FCC Approved, RoHS Compliant keyboard.
- Mechanical sissor switch keys for enhanced tactile feel and reliability.
- Can be supplied without or without Antibacterial silicone coated laser mouse.
- Keyboard Connection USB port with 1.5m cable length.
- Compatible with Microsoft Windows98/Me/2000/NT/XP/Windows 7/8/10

QUICK GUIDE



SterileFLAT Keyboards are quick and easy to clean, with a soft tactile feel. CE, FCC Approved and RoHS compliant.

- **Conventional keyboard design** quickly wipes clean with Antibacterial and Alcohol based agents or Chlor-clean a chlorine based agent as used in many Hospitals.
- **Five levels** of backlit brightness.
- **Low profile Keys** with angled edges and wide caps to aid effective cleaning.
- **Keyboard face Antibacterial**
The membrane combines active Nano Silver technology with a flexible but strong silicone barrier cover.
- **Fully sealed to IP-68 standard.**
- **Mechanical sissor switch keys** for enhanced tactile feel and reliability.
- **Very Cost effective.** Can be supplied with or without antibacterial silicone coated laser mouse.
- **SterileFLAT Keyboards** promote good cross infection practices.
- **Dimensions:** (L x W x H)
385mm x 140mm x 15mm



**Part Number
SF20BL**

Agile Medical

Estate Road 2,
South Humberside Industrial Estate,
Grimsby, DN31 2TG

+44 (0)1472 266 190
sales@agilemedical.co.uk
agilemedical.co.uk

