


Better **solutions.**
Better **care.**

K10 Medical Grade Keyboard



- Complete protection against dust and liquids IP68
- Washable, sterilized in just seconds, no seams or crevices
- Resistant against most common disinfectants and detergents
- CleanRemind (operating time controlled cleaning/ disinfecting reminder)
- Low force and precise feedback key technology
- Rugged shock protected chassis, low noise typing
- EMC tested according EN 60601-1-2
- Compact size, ergonomically low profile design, damage protected key design
- Integrated high resolution capacitive touchpad
- CleanLock-function to avoid unintended inputs
- Magnetic fixing option

Switching technology: Carbon pills on gold contact

Number of Keys: 105 (incl. 12 function keys)

Protection: IP68 (NEMA 4X)

Key travel: 1,25 mm

Actuation force: 1,2 N

Life: approx. 2 Mio. actuations

Chassis: desktop (silicone rubber surface)

Interface: USB

Power: 20mA@5VDC

Dimensions: 390 x 149 x 14 mm

Weight: 1 kg

Temperature: 0°C to +70°C

Humidity: up to 100%

Shock: 3x11 ms pulses (50g) on each of 3 axes

Vibration: 0,04g/Hz, 20Hz-2kHz, 3 hours

Certifications: CE, FCC Class15, Part B, EN 60601-1-2