


One-piece, moulded ABS top  
with fluid retaining edges



HTM71 modular storage trays  
for clear and organised  
management of contents



Smooth surfaces and  
hospital grade materials  
assist cleaning




Tray support system with  
positive stops to prevent trays  
accidentally being withdrawn



Durable rubber buffers  
protect both the trolley and  
fabric of the building




Robust welded steel frame with durable  
epoxy coating maximises product lifetime



**30%** more  
efficient  
storage\*

\*Compared to market alternatives.

**clinell** approved



**Storage Configuration**

5x Medium Polycarbonate Tray - 600W x 400D x 110H (mm)



**Dimensions**

664W x 458D x 950H (mm)

**Side colour options:**



Blue



Red



Green



Yellow




Orange



White


**Tray accessories:**



Identification labels



Handles



Dividers

**Trolley accessories:**



Gel Bottle Holder



Catheter Holder



O<sub>2</sub> Bottle Holder



Glove Box Holder



Sharps Bin Holder



Waste Bin Holder



IV Pole



ISO90012015



NHS approved supplier



CO<sub>2</sub>e Neutral Operations