



STC300

Single Sterile Storage Cart

- Loadable from one side, this cart features a seamless frame construction of stainless steel square tubing with hooks to accommodate 3 to 4 STE baskets or accessories
- Each hook has a weight capacity of 30 kg

Technical description:

- Casters - 4 swivel casters, 2 with double brake, Ø 75 mm, with precision ball bearing
- Tread - Thermoplastic conductive rubber profile, wear-resistant, non-marking
- Caster housing - Zinc-chromated, with bumpers



No. of hook sets	Hook spacing mm	Dimensions			Order code
		Width mm	Depth mm	Height mm	
3	400	610	475	1463	STC300
4	320	610	475	1623	STC400



STC600

Double Sterile Storage Cart

- Loadable from both sides, this cart features a seamless frame construction of stainless steel square tubing with hooks to accommodate 2 x 3 or 2 x 4 sterile goods baskets or accessories
- Each hook has a weight capacity of 30 kg

Technical description:

- Casters - 4 swivel casters, 2 with double brake, Ø 75 mm, with precision ball bearing
- Tread - Thermoplastic conductive rubber profile, wear-resistant, non-marking
- Caster housing - Zinc-chromated, with bumpers



No. of hook sets	Hook spacing mm	Dimensions			Order code
		Width mm	Depth mm	Height mm	
6	400	610	670	1463	STC600
8	320	610	670	1623	STC800