



Better **solutions.**  
Better **care.**

**Agile** Medical  
**Theatre** Essentials

# Agile Medical

Better **solutions.**  
Better **care.**



CO<sub>2</sub>e Neutral  
Operations

Helping to combat climate change  
and sustain our environment for  
future generations.



# Welcome






































Agile Medical is dedicated to improving workflow and efficiency in healthcare environments. Through our wide range of UK manufactured products and a portfolio of digital applications, Agile Medical has been supporting mobilised patient care in the NHS for over 25 years.

Products with the  icon are available to order through NHS Supply Chain eDirect.


SCAN OR CLICK TO VIEW OUR FULL RANGE OF PRODUCTS FOR THEATRES



## In this brochure...

 Theatre Connect	6	 Mayo Tables	28
 Case Study - Manchester NHS Trust	8	 Surgical Steps	29
 Narrow Anaesthetic Trolley	10	 Bowl Stands	30
 LX Anaesthetic Trolley	11	 Sterile Storage Carts	31
 Double Anaesthetic Trolley	12	 Infusion Stands	32
 Carebuddy Trolley	13	 Mobile Theatre Cabinet	33
 ScanCell@	14	 Standard Linen Trolleys	34
 Chart Trolley	14	 Mobile Theatre Trolley	35
 Procedure Trolley	16	 Theatre Stools	36
 Difficult Airway Trolley	17	 Theatre Operator Stools	37
 Infinity Storage	18	 Tech Chairs	38
 Soiled Linen Trolleys	20	 Mobile Privacy Screens	39
 Sackholders	21	 Monitor Arms	40
 Trolley Covers	22	 Heavy Duty Roll Stand	42
 Instrument Trolley with Drawers	23	 Anaesthetic Drug Labels and Tapes	43
 Instrument Trolley - Two Shelves	24	 Pathoban Keyboards and Mice	44
 Instrument Trolley - One Shelf	26	 Theatre Trolley Accessories	45
 Surgical Trolley	27	 Theatre Drawer Accessories	48
		 NHS NHS How to Order	49

Reduce your readiness preparation time by **up to 90%** with Theatre Connect.




**Features:**

- Live trolley readiness
- Automated alerts
- 100% digital process (paperless)
- RFID enabled readiness checks\*
- Early warning system for expiring products
- Automated stock management
- Instantaneous audit data

**Benefits:**

- Theatre trolley checks in under 2 minutes
- Free up medically trained staff time
- Reduce stock wastage
- Quickly identify and locate expired or recalled items
- Real-time compliance dashboard
- Intelligent stock management

**Optimising Theatre Workflow**



- 1** At the end of each day, a digital contents check is performed on all used anaesthetic trolleys using **Theatre Connect**.  
An automated report of missing items for each trolley is sent through to the central store.
- 2** Staff member from the central store then picks all replacement items onto the 'top-up' trolley.
- 3** The 'top-up' trolley is then taken to restock all used anaesthetic trolleys, ensuring 100% readiness in preparation for the next day.

Book a **no-obligation trial** now. Contact [fergus.glass@agilemedical.co.uk](mailto:fergus.glass@agilemedical.co.uk)

# Real Life Operational Benefits:

## A case study of Manchester NHS Trust

“Our Anaesthetic trolley provides the option to anaesthetise in theatre, removing any unnecessary patient movement and improving safety.

The standardised layout and content of our trolleys make it easy for practitioners to transfer between theatres and hospitals and we have eliminated a large amount of waste and out of date products by removing the majority of stock from our Anaesthetic rooms”

**Adam Dobson** *Manchester University NHSFT*



### Patient Safety

The option to anaesthetise in theatre removes unnecessary patient movement to theatre.



### Compliance

Reduce the risk of out of date stock with lower stock levels in fewer locations.



### Practitioner Comfort

Standardised trolley layouts allows practitioners to feel at ease quickly when working in unfamiliar theatres.



### Operational Efficiency

Save up to 15 minutes per patient with equipment more readily available within the operating theatre.



### Stock Control

Reduce stock held in anaesthetic rooms by up to 60% and minimise waste caused by out of date items.



# Narrow Anaesthetic Trolley

Mobilising the storage of anaesthetic equipment reduces staff downtime and prevents 'dead stock' building up. The narrow anaesthetic trolley provides these benefits in a compact, space optimising solution that can be configured to your specific requirements.



### Supplied accessories:



### Trolley Options



Code **NC3405\_ (1 or 2)**  
 Dimensions: 520w, 540d, 946h (mm)  
 5 drawer layout  
 (2x 3", 2x 6", 1x 9")

Colours: **1** White  
**2** Light Blue



Code **NC34GAR\_ (1 or 2)**  
 Dimensions: 520w, 540d, 946h (mm)  
 5 drawer layout  
 (2x 3", 2x 6", 1x 9")

Colours:  White Paediatric (Garden Theme)  
 Light Blue Paediatric (Garden Theme)

For all accessories see pages 45-48.

# LX Anaesthetic Trolley

Tried and tested in the theatre departments of many leading NHS hospitals, this range exceeds the infection control and durability requirements of a theatre environment.



### Trolley Options



Code **LX31ANE\_ (1 or 2)**  
 Dimensions: 690w, 540d, 930h (mm)  
 4 drawer layout  
 (1x 3", 2x 6", 1x 9")

Colours: **1** White  
**2** Light Blue



Code **LX34ANE\_ (1, 2 or 3)**  
 Dimensions: 690w, 540d, 1010h (mm)  
 5 drawer layout  
 (2x 3", 2x 6", 1x 9")

Colours: **1** White  
**2** Light Blue  
**3** Blue




Code **LX34GAR\_ (1 or 2)**  
 Dimensions: 690w, 540d, 1010h (mm)  
 5 drawer layout  
 (2x 3", 2x 6", 1x 9")

Colours:  White Paediatric (Garden Theme)  
 Light Blue Paediatric (Garden Theme)

For all accessories see pages 45-48.

### Supplied accessories:



# Double Anaesthetic Trolley

The optimum storage solution for all your anaesthetic equipment. The combination of LX and NC storage on one trolley brings greater flexibility with equipment storage and ensures all the vital equipment is right at hand throughout the entirety of the procedure.



**Supplied accessories:**

- Slide pull-out extension table
- 10 x Label Holders
- 4 x Divider sets
- Push handle



**Trolley Options**



Code **DC34092**  
Dimensions: 1125w, 545d, 1010h (mm)

Colours:  Light Blue



Code **DC34GAR\_ (1 or 2)**  
Dimensions: 1125w, 545d, 1010h (mm)

- Colours:
-  White Paediatric (Garden Theme)
  -  Light Blue Paediatric (Garden Theme)

VIEW OUR FULL RANGE OF PRODUCTS FOR THEATRES



For all accessories see pages 45-48.

# Carebuddy Trolley

Designed and manufactured in the UK, the Carebuddy range of trolleys from Agile Medical provides the optimum solution for mobilising storage of medical supplies throughout healthcare facilities.



Various sizes, colours and configurations possible (contact us today to configure your Carebuddy).

To order add colour code required (e.g; CBOB10BL050)



**Trolley Options**



Code **CBOB10**  
Dimensions: 664w, 458d, 950h (mm)

**Storage configuration:**  
5x Medium Polycarbonate Tray - 600w, 400d, 110h (mm)

**Tray accessories:**



Identification labels



Handles



Dividers

**Trolley accessories:**

- Gel bottle holder
- Catheter holder
- O<sub>2</sub> bottle holder
- Glove box holder
- Sharps bin holder
- Waste bin holder
- IV pole



**30%**  
more efficient storage\*

\*Compared to market alternatives

# ScanCell®

The high quality mobile storage solution.



**ScanCell® is our range of high quality mobile storage solutions that utilises our modular system of trays and baskets.**

### Key benefits

#### Flexible

Highly customisable to meet specific job functions and needs

#### Light weight

Light weight trolleys for easy handling

#### Ergonomic

Job specific trolleys ensure an improved work environment

#### Easy

High quality swivel wheels for increased maneuverability.

**CLICK OR SCAN THE QR CODE TO CONFIGURE YOUR SCANCELL® TROLLEY NOW**



# Chart Trolley

The addition of an integrated chart board to our most popular storage trolley range provides the optimum solution for maintaining patient records in recovery areas. The chart board can be adjusted to accommodate all users and includes two bulldog clips for securing the charts and a tray for storing pens.



### Supplied accessories:

Slide pull-out extension table



4 x Label Holders



2 x Divider sets



Push handle



Chart board



### Trolley Options



Code **NC3405\_/CT (1 or 2)**  
Dimensions: 520w, 540d, 946h (mm)

**Colours:** **1** White  
**2** Light Blue



Code **LXC3405\_ (1 or 2)**  
Dimensions: 690w, 540d, 1010h (mm)

**Colours:** **1** White  
**2** Light Blue



Code **DC3409\_/CT (1 or 2)**  
Dimensions: 1125w, 545d, 1010h (mm)

**Colours:** **1** White  
**2** Light Blue

**VIEW OUR FULL RANGE OF PRODUCTS FOR THEATRES**



**SCAN OR CLICK ME**

For all accessories see pages 45-48.



# Procedure Trolley

Our Procedure Trolleys can be arranged to suit a multitude of medical departments and situations.



### Supplied accessories:

- Slide pull-out extension table
- 4 x Label Holders
- 2 x Divider sets
- Push handle



### Trolley Options



Code **LX28PRO\_ (1 or 2)**  
 Dimensions: 690w, 540d, 850h (mm)  
 4 drawer layout  
 (1x 3", 3x 6")  
**Fits under standard worktop**

- Colours:
- 1** White
  - 2** Light Blue



Code **LX31PRO\_ (1, 2 or 3)**  
 Dimensions: 690w, 540d, 930h (mm)  
 4 drawer layout  
 (1x 3", 2x 6", 1x 9")

- Colours:
- 1** White
  - 2** Light Blue



Code **LX34PRO\_ (1 or 2)**  
 Dimensions: 690w, 540d, 1010h (mm)  
 5 drawer layout  
 (2x 3", 2x 6", 1x 9")

- Colours:
- 1** White
  - 2** Light Blue

For all accessories see pages 45-48.

# Difficult Airway Trolley

The Difficult Airway Trolley can be customised, with 4 or 5 drawer option available and a zero risk emergency breakaway locking system and many other unique features. Each trolley is supplied with drawer front labels. These labels help clinical and non-clinical staff identify the correct procedure in a difficult airway or intubation situation.



### Supplied accessories:

- Breakaway Seals
- 5 x Label Holders
- Slide pull-out extension table



### Trolley Options



Code **LX28AIR\_ (5 or 8)**  
 Dimensions: 690w, 540d, 850h (mm)  
 4 drawer layout  
 (1x 3", 3x 6")  
**Fits under standard worktop**

- Colours:
- 5** Red
  - 8** Dark Green



Code **LX31AIR\_ (5 or 8)**  
 Dimensions: 690w, 540d, 930h (mm)  
 4 drawer layout  
 (4x 6")

- Colours:
- 5** Red
  - 8** Dark Green



Code **LX34AIR\_ (5 or 8)**  
 Dimensions: 690w, 540d, 1010h (mm)  
 4 drawer layout  
 (3x 6", 1x 9")

- Colours:
- 5** Red
  - 8** Dark Green

For all accessories see pages 45-48.

# Infinity Storage

## HTM71 Materials Management Modular Storage.

Agile Medical's innovative, high-quality aluminium product is a flexible system that expands easily to different widths and depths. Combined with a wide range of accessories, this is a versatile storage system that satisfies all requirements.

Create your ideal **storage and logistics environment**



From shelf to all-rounder **Storage, preparation, work**

**100% compatible**

Infinite possibilities - **A personalised storage system**



**CLICK OR SCAN THE QR CODE TO CONFIGURE YOUR INFINITY STORAGE SYSTEM NOW**



**Ideal working conditions**  
Keep your products clearly sorted, stored safely, and within easy reach. This saves both time and space! Coupled with our comprehensive labelling system, you can minimise errors while boosting efficiency.

# Soiled Linen Trolleys

- Lightweight, anodised aluminium body
- One-piece, moulded ABS base and lid
- Bag-holder with front opening for easy removal of linen bags
- Bag grip-ring with rubber corners to hold linen bags in place
- 4 swivel castors (Ø75mm)
- Integrated foot operated lid mechanism
- Suitable for disinfection
- (Linen sacks not included)



**Lid colour options:**

White, Green, Blue, Red, Yellow, Brown, Grey, Black



## Trolley Options



**Single Soiled Linen Trolley**  
Code **STS348**

Dimensions: 390w, 410d, 880h (mm)

Lids to be ordered separately.



**Double Soiled Linen Trolley**  
Code **DTD748**

Dimensions: 785w, 410d, 880h (mm)

Lids to be ordered separately.



**Triple Soiled Linen Trolley**  
Code **TTT114**

Dimensions: 1180w, 410d, 880h (mm)

Lids to be ordered separately.

# Sackholders

- 100% durable plastic construction
- Moulded handsfree frame with foot pedal operation
- Removable body
- Easy clean surfaces and radius corners
- 100% corrosion resistant foot unit
- Available in all standard waste segregation colours



**Lid colour options:**

White, Yellow, Orange, Black, Tiger Stripe



**28ltr Capacity**  
Code **AM28LWHWH**

Dimensions: 250w, 355d, 675h (mm)



**42ltr Capacity**  
Code **AM42LWHWH**

Dimensions: 485w, 390d, 570h (mm)



**80ltr Capacity**  
Code **AM80LWHWH**

Dimensions: 485w, 390d, 847h (mm)

**VIEW OUR  
FULL RANGE  
OF PRODUCTS  
FOR THEATRES**



**SCAN OR  
CLICK ME**

## Trolley Covers

In a busy hospital environment, it can be difficult to keep your trolleys and equipment clean when not in use. That is why Agile Medical has developed a wide range of practical Trolley and Equipment Covers.

Our Trolley and Equipment Covers are manufactured from high quality, wipe-clean PVC and feature a choice of premium zip and/or velcro fastenings.

Our design team engaged in rigorous testing for this product, creating a practical front edge fastening to simplify the removal process and ensure your weekly equipment clean is a quick and hassle-free experience.



Our transparent Trolley & Equipment Covers allow your equipment to remain easily identifiable whilst in storage, while simultaneously preventing the build-up of dust, dirt and infectious particles that may harm your patients.

### Features

- Clear PVC design
- Large range of sizes available
- Custom sizing available
- Choice of zip and/or velcro fastenings
- Wipe clean.



### Custom Design Service Available



## Instrument Trolley with Drawers

Slim-line Instrument Trolley with 2 storage drawers and 2 storage shelves, ideal for hospitals and health care centres. Optional moulded ABS top with fluid retaining edges.



### Trolley Options



Code **IT3502SW**  
Dimensions: 500w, 500d, 900h (mm)



**Single ABS Top**  
Code **IT5050T**  
Dimensions: 500w, 500d (mm)

VIEW OUR  
FULL RANGE  
OF PRODUCTS  
FOR THEATRES



SCAN OR  
CLICK ME

# Instrument Trolley - Two Shelves

Instrument Trolley – Strong (grade 304) stainless steel seamless welded construction with rounded corners (no sharp edges), swivel anti-static castors and rubber buffers.



**Features:**

- Anti-static castors
- Flush-welded, smooth finish and easy clean
- Buffers



**Trolley Options**



Code **SS4545F**  
Dimensions: 450w, 450d, 840h (mm)



Code **SS7545F**  
Dimensions: 750w, 450d, 840h (mm)



Code **SS9145F**  
Dimensions: 910w, 450d, 840h (mm)

Instrument Trolley – Strong (grade 304) stainless steel seamless welded construction. Flush welded shelves with raised rim containment edges, swivel anti-static castors and rubber buffers.



**Features:**

- Anti-static castors
- Flush welded shelves
- Buffers



**Trolley Options**



Code **SS4545R**  
Dimensions: 450w, 450d, 840h (mm)



Code **SS6145R**  
Dimensions: 610w, 450d, 840h (mm)



Code **SS9145R**  
Dimensions: 910w, 450d, 840h (mm)

## Instrument Trolley - One Shelf

Tall Instrument Trolley – Strong (grade 304) stainless steel seamless welded construction with rounded corners (no sharp edges), swivel / lockable castors. Available for 900mm and 1000mm workshop heights.



### Trolley Options



Code **SIT0994**  
Dimensions: 900w, 450d, 900h (mm)



Code **SIT1094**  
Dimensions: 900w, 450d, 1000h (mm)

## Surgical Trolley

Special Theatre Trolley – Strong (grade 304) stainless steel seamless welded construction. All shelves with raised rim containment edges, swivel / lockable castors and removable 15L stainless steel buckets (available separately).



### Options



**Surgical Trolley**  
Code **CIT0947**  
Dimensions: 770w, 450d, 900h (mm)



**15L Bucket**  
Code **NSB300**  
Diameter: 295mm  
Height: 310mm

VIEW OUR  
FULL RANGE  
OF PRODUCTS  
FOR THEATRES



SCAN OR  
CLICK ME

# Mayo Tables

Easy to use, mobile Mayo Tables providing a place for sterile instruments and supplies used during surgery. Fully stainless-steel welded construction providing compatibility with all standard hospital cleaning products.



## Options



**Mayo Table with manual height adjustment**  
Code **HMTD46**

- Stainless steel Frame (AISI 304)
- ST/ST top with perimetrical rim (50cm x 70cm WxD)
- Manual height adjustment (80-120 cm)
- 4 medical grade castors (2 braked)



**Mayo Table with hydraulic height adjustment**  
Code **HMTD48**

- Stainless steel Frame (AISI 304)
- Removable ST/ST top with perimetrical rim (50cm x 70cm WxD)
- Hydraulic foot pedal height adjustment (80-120 cm)
- 4 medical grade castors (2 braked)



**Mayo Table with hydraulic height adjustment**  
Code **HMTD50**

- Stainless steel Frame (AISI 304)
- ST/ST top with perimetrical rim (50cm x 70cm WxD)
- Hydraulic foot pedal height adjustment (80-120 cm)
- 4 medical grade castors (2 braked)

# Surgical Steps

Surgeon Steps – Stable 4-footed single step construction, grade 304 stainless steel seamless welded frame construction with removable stainless steel ribbed tread for easy cleaning. 2 non-slip, non-marking adjustable rubber feet and 2 conductive polymer screw-in feet.



## Options



Code **STE1631**  
Dimensions: 600w, 350d, 130h (mm)



Code **STE1632**  
Dimensions: 600w, 350d, 220h (mm)



Code **STE1633**  
Dimensions: 600w, 350d, 300h (mm)



Code **STE26624**  
Dimensions: 600w, 650d, 220/440h (mm)

## Features:

2 x Non-slip Adjustable Rubber Feet



Removable Top (S/Steel Non-Slip)



2 x Conductive Polymer Screw-in Feet



Ribbed Tread for Easy Cleaning



## Bowl Stands

Manufactured from grade 304 stainless steel, our bowl stands and kickabouts exceed infection control requirements for within sterile environments. All parts can easily be cleaned and sanitised.



### Options



**Single Bowl Stand**  
Code **AM/BSD-79**

Dimensions: 350w, 350d, 800h (mm)

- Stainless steel Frame (AISI 304)
- 5-star ST/ST base with swivel castors (2 braked)
- 1 removable ST/ST bowl



**Double Bowl Stand**  
Code **AM/BSD-80**

Dimensions: 700w, 350d, 800h (mm)

- Stainless steel Frame (AISI 304)
- 5-star ST/ST base with swivel castors (2 braked)
- 2 removable ST/ST bowls

## Sterile Storage Carts

Single & Double Sided Storage Carts – Strong (grade 304) stainless steel seamless welded frame with hook configuration to accommodate 3 – 4 baskets, hooks tested to 50kg load capacity, swivel / lockable castors.



### Options



**Single Storage Cart**  
Code **STC300**

Dimensions: 610w, 470d, 1450h (mm)



**Double Storage Cart**  
Code **STC600**

Dimensions: 610w, 680d, 1450h (mm)

VIEW OUR  
FULL RANGE  
OF PRODUCTS  
FOR THEATRES



SCAN OR  
CLICK ME



# Infusion Stands

Light, medium and heavy duty Infusion Stands to suit a variety of applications – Robust design, stable 5-star base, easy to clean and simple to operate.



## Options



### Light Duty Infusion Stand

2 Hooks **NS17230**  
4 Hooks **NS17229**

- Colour of Plastic Parts: Blue
- Material: Stainless Steel
- Height Adjustment: 1350-2150mm
- Base Diameter: 635mm
- Base Weight: 3.5kg



### Medium Duty Infusion Stand

2 Hooks **ICS17230**  
4 Hooks **ICS17229**

- Colour of Plastic Parts: Blue
- Material: Stainless Steel
- Height Adjustment: 1500-2350mm
- Base Diameter: 680mm
- Base Weight: 8.5kg / 6.5kg



### Heavy Duty Infusion Stand

Code **PCI172954**

- Specifically for pump attachments
- 4 Hooks – Stainless Steel
- Colour of Plastic Parts: Black
- 5 star 630mm diameter base
- Height Adjustment: 1430-2260mm (soft drop system)
- Load capacity: 32kg
- Twin castors - electrically conductive and lockable

# Mobile Theatre Cabinet

Single or Double HTM71 Mobile Base Cabinets are supplied with lockable doors and internals. Specially designed to fit below 900mm worktop heights.



## Cabinet Options



### Single B40 Cabinet

Code **SMC846\_**

Dimensions: 450w, 650d, 860h (mm)

- Comes with:
- 1 x Door (lockable)
  - 4 x 600x400x100mm Poly Trays
  - 1 x 600x400x200mm Poly Basket



### Double B40 Cabinet


Code **DMC896\_**

Dimensions: 900w, 650d, 860h (mm)

- Comes with:
- 2 x Doors (lockable)
  - 8 x 600x400x100mm Poly Trays
  - 2 x 600x400x200mm Poly Basket

The above cabinets are suited for 600mm insertions. Also available to suit 400mm insertions.

Type A - 600mm insertion



Type B - 400mm insertion



### Various door colours available (contact us for full list):

To order add colour code required (e.g; DMC896-**U4817**)



# Standard Linen Trolleys

A range of easy-to-clean, mobile trolleys to provide enclosed storage of clean linen. One of our best selling ranges for Endoscopy, Theatres and Day Surgery.



## Trolley Options



**TYPE A**  
Code **LT1085**  
Dimensions: 740w, 535d, 1000h (mm)  
Comes with 1 central shelf.



**TYPE B**  
Code **LT1285**  
Dimensions: 740w, 535d, 1190h (mm)  
Comes with 2 shelves.



**TYPE C**  
Code **LT1485**  
Dimensions: 740w, 535d, 1400h (mm)  
Comes with 2 shelves.

## Double Linen Cabinet

Code **LTC875**  
Dimensions: 700w, 540d, 800h (mm)

Comes with two internal compartments with adjustable shelves.



# Mobile Theatre Trolley

An attractive and specially designed HTM71 mobile cabinet for the theatre departments. High impact strength, easy-to-clean and resistant to chemicals and high energy radiation (gamma and x-rays).



## Trolley Options



**Mobile Theatre Trolley**  
Code **MTC1675\_**  
Dimensions: 740w, 580d, 1655h (mm)  
Body: Stainless Steel  
Door Frame: PE500 with part transparent  
Capacity:  
- 9 x 600x400x100mm Trays  
or  
- 8 x 600x400x100mm Trays  
- 1 x 600x400x200mm Tray

VIEW OUR FULL RANGE OF PRODUCTS FOR THEATRES



SCAN OR CLICK ME

## Various door colours available:

To order add colour code required (e.g; MTC1675-6)



# Theatre Stools

Medical Line range with comfortable 40cm round seat, attractive colour options and 10cm height adjustment creating an optimum working area of between 80-90cm.



## Options



**Medical Stool (Round Seat)**  
Code **MR6101\_**  
Height adjustment 570-760mm  
Seat diameter 400mm



**Medical Stool (Round Seat) with backrest**  
Fixed seat angle **MR6211F\_**  
Adjustable seat angle **MR6211T\_**  
Height adjustment 570-760mm  
Seat diameter 400mm



**Medical Stool (Round Seat) with backrest**  
Fixed seat angle **MR6221F\_**  
Adjustable seat angle **MR6221T\_**  
Height adjustment 570-760mm  
Seat diameter 400mm



**Height Adjustment by Foot**  
Code **FA6000**

All models are also available with the **Ergo seat**. Contact us for product details.



### Various colour fabrics available.

To order add colour code required (e.g; MR6211F-**K89**)



# Theatre Operator Stools

Our range of operator stools designed and developed with input from leading ergonomists and surgeons to support theatre operators throughout procedures, by creating a stable and ergonomic sitting posture.



## Options



**Amazon Saddle Stool**  
Fixed seat angle **MAS6311F\_**  
Adjustable seat angle **MAS6311T\_**  
Height adj. 590-770mm  
Seat appr. 340w, 300d (mm)



**Jumper Saddle Stool**  
Fixed seat angle **MJS6311F\_**  
Adjustable seat angle **MJS6311T\_**  
Height adj. 590-770mm  
Seat appr. 440w, 300d (mm)



**Ergo Support Surgeons' Stool**  
Code **ME6331TFHB**  
- Adjustable seat angle  
- Foot operated height adj. 520-710mm  
- Central foot brake  
- 2 swivel, adjustable armrests



**Height Adjustment by Foot**  
Code **FA6000**

### Various colour fabrics available.

To order add colour code required (e.g; MAS6311-**K89**)



## Tech Chairs

- PU Seat and Back Upholstery (easy to clean), Dark Grey
- Hydraulic Hand Operated Height Adjustment
- Independent Back Height and Back Angle Adjustment
- 60cm Synthetic Base, Black
- 5 Soft Swivel Castors
- Optional Armrests



### Options



**Option 1**  
Code **SPR1500**

No footrest  
44-57cm Adjustment



**Option 2**  
Code **SPR1501**

With footrest  
52-71cm Adjustment

## Mobile Privacy Screens

Quick and easy to open out and close away, these mobile privacy screens are ideal for having ready-to-hand in a multi-bay recovery area to help maintain the privacy of patients.



### Options



**3 Panel Privacy Screen**  
Code **AM31400**

Dimensions: 1650L x 1680H (mm)


- 3 white multi-use vinyl panels
- Aluminium tube frame
- Swivel Castors Ø50mm



**4 Panel Privacy Screen**  
Code **AM31401**

Dimensions: 2220L x 1680H (mm)


- 4 white multi-use vinyl panels
- Aluminium tube frame
- Swivel Castors Ø50mm



# Monitor Arms

## Mounting solutions that set new standards.

Compliant with the highest medical regulations, plus integrated cable routing and a 5-year maintenance-free warranty. What's really at our core is Infection Prevention. CIM arms are clinically proven to be the best in the world for Infection Prevention, according to independent test results.



### Infection prevention through easy-to-clean surfaces.

Clinical test results, comparing all of the leading brands in medical mounting solutions, revealed CIM arms to be at least four times better for advanced hygiene control.



### Integrated cable management

We meet the growing demands of hygiene and aesthetics in the medical environment while protecting cables from damage and accidental disconnection.



### Silky ergonomics

Smooth adjustment and maneuverability with perfect placement of the monitors and devices.



### Flexible positioning

Large swiveling and turning ranges for optimum positioning of the devices.



### Space saving

Optimal parking positions through flexibility.



### Zero maintenance

Guaranteed maintenance-free operation during the standard 5-year warranty.



### 6-fold safety factor

The highest safety standard on the market for patient, user and device safety.



### TÜV-certified safety

Independently tested systems that comply with current EN 60601 regulations.



### Precise engineering

Guarantees a level screen in all viewing positions. Finesse at your fingertips.



**CIM mounting solutions are compatible with patient monitors from a wide range of medical device manufacturers.**

To discuss the various versions, colour schemes, and mounting options available, please get in touch with our product support team at [sales@agilemedical.co.uk](mailto:sales@agilemedical.co.uk).

# Heavy Duty Roll Stand

The Heavy-Duty Roll Stand can be configured according to your needs with an extensive range of accessories available to match your requirements. All individual solutions are tested and certified before delivery and are tip tested in conformance with EN 60601-1, 3rd Edition.



**Accepts:**

**Drawers and shelves:** Single and double drawer consoles or simply shelves to place your equipment.

**Support arms:** CIM height adjustable and non-height adjustable arms to mount monitors, devices and other medical equipment.

**Accessories:** A huge variety of accessories is available: cable hooks, baskets, IV poles.

## Monitor Presentations

Which monitor would you like to mount? Agile Medical offers solutions for a wide range of medical device manufacturers.

<p><b>A</b></p>  <p>5" Plunge plate</p>	<p><b>PH</b></p>  <p>Philips</p>
<p><b>MD</b></p>  <p>Mindray</p>	<p><b>FK</b></p>  <p>Fukuda Denshi</p>
<p><b>SLX</b></p>  <p>Spacelabs Xprezzon</p>	<p><b>G</b></p>  <p>VESA 75/10</p>
<p><b>Q</b></p>  <p>VESA quick release</p>	<p><b>TN</b></p>  <p>Laptop</p>

# Anaesthetic Drug Labels and Tapes

Helping caregivers correctly communicate injectable medications to prevent confusion or mistakes in administration.



**Meet critical Patient Safety & Regulatory Guidelines with accurate medication labelling.**

Agile Medical's Anaesthesia Labels & Tapes help prevent medication errors caused by the use of unlabelled or mislabelled medications and other procedural breakdowns that can have a fatal impact on patients. Safe medication practices are critical to patient safety and preventing adverse drug events.

- Accurate medication identification
- Compliance with specific international guidelines and regulatory requirements
- ASTM standardised colours, drug classifications and required identification fields

**Contact us today to receive your drug label order form.**



**Need tailored? We can print customised labels.**

# Pathoban Keyboards and Mice



**Pathoban USB Wired Keyboard**  
Code **SK10801**

**Pathoban Wireless Medical Silicone Keyboard**  
Code **SK108005**

- Features:**
- IP classification: IP68 water and dust protection
  - Language layouts available on request
  - With shortcut key of Calculator, Print, Lock keyboard, Cleaning mode
  - Wired version with waterproof Silicone USB cap
  - 100% Silicone body (Latex free)
  - Real scissor type keys: with excellent tactile responsiveness
  - Plug in and play. No drivers needed
  - With cleaning lock key and Cleaning reminder function
  - Only 0.5mm raised surfaces: easy-to-clean and keys positioned
  - Suitable for all kinds of abluent and bactericides.

**Pathoban Sealed Washable Mouse**  
Code **ASM11**

- Features:**
- Fully sealed
  - IP68 water and dust protection, chemical resistant
  - Durable silicone rubber material
  - Capacitive touch sensor in middle area, 3 standard mouse left key and right key, middle key
  - Support to customise logo
  - Color available: White/Black
  - Wired and wireless versions are all available.



# Theatre Trolley Accessories



**Anaesthesia Label Dispenser**  
Clear plastic anesthesia label dispenser.  
Dimensions: 300 x 127 x 140mm  
Code: **AMDSP-HC10**



**Horizontal 3 box glove holder – clear ABS**  
Suitable for narrow trolleys.  
Dimensions: 260 x 395 x 105mm  
Code: **GH2310W**



**Vertical 3 box glove holder – clear ABS**  
Dimensions: 400 x 258 x 105mm  
Code: **GH4210W**



**IV pole**  
Adjustable IV pole  
Code: **IVP100**



**Tilt Bin Organiser (Half Size)**  
Suitable for narrow trolleys.  
Dimensions: 300 x 111 x 230mm  
Code: **TBS312**



**Tilt bin organiser (small)**  
600 x 111 x 230mm  
Code: **TBS612**



**Tilt bin organiser (large)**  
600 x 161 x 168mm  
Code: **TBL611**



**Tilt bin organiser (extra large)**  
600 x 200 x 168mm  
Code: **TBX612**



**Accessory bridge**  
800mm  
Code: **AB800**



**Glove dispenser & back rail**  
160 x 260 x 95mm  
Code: **GDB119**



**3 box glove holder & side rail**  
400 x 260 x 85 xmm  
Code: **GDS428**



**Bulk storage bins (L) & back rail**  
75 x 105 x 190mm  
Code: **BSB711**



**Glove dispenser & side rail**  
160 x 260 x 95mm  
Code: **GDS119**

# Theatre Trolley Accessories



**Waste bin & side rail**  
265 x 255 x 460mm  
Code: **WBS224**



**Waste container w/lid & side rail**  
255 x 170 x 315mm  
Code: **WCS213**



**Garbage bin & side rail**  
255 x 170 x 315mm  
Code: **GBS213**



**Oxygen holder (L) & side rail**  
125 x 125 x 700mm  
(E Size)  
Code: **OTL117**



**Suction Machine Holder & side rail**  
80 x 400 x 200mm  
Code: **SMH300**



**Catheter holder & side rail**  
100 x 100 x 480mm  
Code: **CHS114**



**Multi-storage & side rail**  
Code: **MSL100**



**File holder & side rail**  
571 x 397 x 150mm  
Code: **FHS531**



**Sharps bin holder & side rail**  
Code: **SBH100**



**Catheter holder B & side rail**  
Type B (Vertical Storage)  
600 x 390 x 150mm (HWD)  
Code: **301106/SR**



**Apron roll dispenser & side rail**  
Can be trolley mounted  
162 x 162 x 410mm  
Code: **306171/SR**



**Monitor arm and laptop holder**  
Code: **LAP100**



**Sharps bin holder & side rail**  
100 x 190 x 110mm  
Code: **SBH200**



**Sharps bin holder**  
350 x 170mm  
Code: **SBH300**



**Gel holder**  
Code: **GMH100**



**Monitor arm and screen holder**  
Code: **COM100**



**Arm support holder**  
Support the laptop or computer arm and fix to the trolley  
Code: **ASH100**



**PC holder and 2 x side rails**  
Code: **PCH400**



**Foldaway side drawer unit**  
ABS 3 side files (open or closed)  
Code: **FAS416**



**Stainless steel defibrillator shelf (extra strength)**  
330 - 410 x 330mm  
Code: **DSS433**



**Oxygen holder (S) & side rail**  
145 x 125 x 403mm  
(CD Size)  
Code: **OTS114**



**Side rails**  
450mm  
Code: **ASR450**



**Back rails**  
LX Trolleys - 580mm  
Code: **ABR550**  
  
DC Double Trolleys - 980mm  
Code: **ABR980**



**Trolley keys**  
Code: **LXK200**



**Magnetic label holders**  
80 x 50 x T5mm  
6pcs/pack  
Code: **MLH855**



# Theatre Drawer Accessories



**Straight divider set**  
76mm fixed divider set in drawer (3 x short + 3 x long + 2 x strip)  
Code: **LXF076**



**Straight divider set**  
152mm fixed divider set in drawer (2 x short + 2 x long + 2 x strip)  
Code: **LXF152**



**Adjustable divider set**  
76mm adjustable divider set in drawer (3 x short + 6 x adjustable + 2 x 76mm strip + 6 x T clips)  
Code: **LXA076**



**Adjustable divider set**  
152mm adjustable divider set in drawer (3 x short + 4 x adjustable + 2 x 152mm strip + 4 x T clips)  
Code: **LXA152**



**Laptop/keybaord drawer (76mm)**  
All colours. Available for all LX trolleys (690mm wide)  
Code: **PCB690**



**ABS trays - small**  
Size: 71 x 570 x 396mm  
Code: **LXT764**




**ABS trays - large**  
Size: 152 x 570 x 396mm  
Code: **LXT164**



**Straight divider set (for trays)**  
76mm fixed divider set (3 x short + 3 x long)  
Code: **FTD076**



**Straight divider set (for trays)**  
152mm fixed divider set (2 x short + 2 x long)  
Code: **FTD076**



**Divider set for narrow trolleys**  
Size: 71 x 399 x 398 mm  
Code: **NAA076**



**Divider set for narrow trolleys**  
Size: 152 x 399 x 398 mm  
Code: **NAA152**



How to Order: 

If the NHS Trust you work at has Agile Medical set up as a supplier, you will just need a quotation from Agile Medical for the products you require.

Depending on how the Trust internal order process works, the quotation should be used to raise a requisition internally for Procurement to raise an official Purchase Order.

 Orders can be sent to our [sales@agilemedical.co.uk](mailto:sales@agilemedical.co.uk) email address.

 For orders over the phone, our representatives will require your full contact details, invoice and delivery address as well as the Purchase order number. A proforma invoice with our company bank details will be sent to your Finance team to make payment.

 Products with the **e** icon are available to order through NHS Supply Chain eDirect.

VIEW OUR FULL RANGE OF PRODUCTS FOR THEATRES





Better **solutions.**  
Better **care.**



ISO9001:2015



NHS approved supplier



CO<sub>2</sub>e Neutral  
Operations