



Better **solutions.**  
Better **care.**

# Agile Medical

## Vena<sup>®</sup> Vario Phlebotomy Chair

Excellent design for the comfort of patient and nurse.

**The Vena® Vario Phlebotomy chair combines comfort and seamless functionality. Easy operation makes it possible to act quickly in emergency situations, with the patient able to be moved safely and quickly into the Trendelenburg position.**

The Vena® Vario phlebotomy chair makes it possible to sit comfortably and relax during blood sampling. The nurse can provide back and leg support simultaneously which ensures good support to the patient's whole body. When reclined, the nurse or medical expert can also use it as a treatment table. It is equipped with two height adjustable armrests for the left and right arm. By opting for the multi-adjustable armrests it is even easier to put the patient's arm in a correct and relaxed position and the nurse can draw blood from the patient's arm in an ergonomically responsible manner.



**Standard model:**

Armrests: 2 x fixed Pur  
 Fixed height: 46cm  
 Sizes: Reclined position 167.5cm. Base 58cm(W) x 70cm(L)  
 Robust frame (white RAL 9016)  
 Protective glides

**Options/accessories:**

Armrests: Multi adjustable in PU or upholstered in Stamskin or PUxx  
 Get-up brackets  
 Hospital grade swivel castors (with brake)

**Stamskin colours**



**PUxx colours**



**Stamskin®** 100% waterproof; Urine-, blood-, sweat-, fungal proof, Clean with water and soap, rinse of with water,

**PUxx®** 100% waterproof, PVC and Phthalate free; resistant against commonly used cleaning products and disinfectants (alcohol based disinfectants, active chlorine, quaternary ammonium compounds).