



Better **solutions.**
Better **care.**

Agile Medical

Vena Vario Levant Phlebotomy Chair

The #1 choice phlebotomy chair for comfortable seating and healthy working

With the Vena® Vario phlebotomy chair, comfort and functionality are seamlessly combined. Easy operation makes it possible to act quickly in emergency situations. The patient can be quickly placed into Trendelenburg-position in just one simple movement.

The design of the Vena Vario ensures a patient can sit comfortably and relaxed during blood sampling. The operator or nurse can adjust the back and leg support synchronously, making it possible to draw blood while a patient is either seated or laid down. The chair is not only comfortable for the patient, but also the operator, who can do his or her job in an ergonomic and healthy manner. The Vena Vario phlebotomy chair is standard equipped with two height adjustable phlebotomy armrests. Optionally, these armrests can be upgraded to multi-adjustable armrests that make it easier to bring the vein of the patient in line with the work position of the nurse. The Vena Vario Levant model is height adjustable with a high class hydraulic foot pump.

Vena Vario Levant standard:

- Armrests: 2 x 'Vena® fixed' with PU arm pads
- Seating height: 47 - 67 cm
- Measurements: reclined position 167.5 cm; base 58 cm (W) x 70 cm (L)
- Robust frame (white RAL 9016)
- Hospital grade castors (with brakes)
- Hand-controlled synchronic back and leg rest operation
- Trendenbureg position 7°

Options and accessories:

- Armrests: 2 x 'Vena® multi adjustable' with PU arm pads
- Armrests: 2 x 'Vena® multi adjustable' with pads upholstered in Stamskin or PUxx
- Get-up brackets
- Central foot brake (on castors front side)
- Foot-control for foot back and leg rest adjustment and Trendelenburg release




Also available: Vena Vario Levant on base

- Fully rotating Vena Vario Levant on a base
- Height adjustment by means of a hydraulic foot pump

Stamskin colours



PUxx colours



Stamskin® 100% waterproof; Urine-, blood-, sweat-, fungal proof, Clean with water and soap, rinse of with water,

PUxx® 100% waterproof, PVC and Phthalate free; resistant against commonly used cleaning products and disinfectants (alcohol based disinfectants, active chlorine, quaternary ammonium compounds).