



Examination/Treatment Furniture

CONTENTS

Medical Couches

Two Section Couches	114
Three Section Couches	116
Paediatric Couches	118
Bariatric Couches	120

Gynaecology Couches

Obstetrics/Gynaecology Couches	122
--------------------------------	-----

Procedure Chairs

Obstetrics/Gynaecology Chairs	124
Multi-Purpose Clinic Chairs	126
Paediatric Multi-Purpose Clinic Chairs	128
Phlebotomy Chairs	130
Vena Vario Phlebotomy Chairs	132
Ultrasound Scanning/ Medical Procedure Chairs	134

Mammography Chair

136

Two Section Couches

Multi-discipline two-section models suit a variety of medical, clinical and minor surgical applications. Electronic variable height comes as standard, with an optional hydraulic lift alternative available on selected couches.



Attachments to order:



*CC2MD/25 & CC2MD/28 only

Features:

- Gas assisted adjustable backrest section
- Retractable wheelbase with foot lever control
- Anti-microbial seamless 2" vinyl upholstery
- Steel frame with white anti-microbial protective coating.

Dimensions:

Length: 1900mm Height adjustment: 480-890mm Safe working load: 225kg

Vinyl Colours:

- | | | | | |
|-------------|-----------|------------|-----------|--------------|
| 1 Burnt Red | 2 Rose | 3 Mulberry | 4 Lilac | 5 Navy |
| 6 Royal | 7 Lido | 8 Aqua | 9 Ocean | 10 Citrus |
| 11 Sage | 12 Grey | 13 Dove | 14 Ginger | 15 Buttercup |
| 16 Taupe | 17 Barley | 18 Stone | 19 White | 20 Black |

Options to order:



Two Section Couches

Ordering:

Add the chosen colour code to the end of the product code.
For example: **CC2MD/25H5** for Rose Two Section Couch - Hydraulic

Product	Code
Two Section Couch - Electronic	CC2MD/25E
Two Section Couch - Hydraulic	CC2MD/25H
Attachments to order	
Hygiene roll holder	HRH001
Safety side rails (pair)	SSR001
Safety side pads (pair)	SSP001
Head cushion support	HCS001
IV pole	IVP001
Ankle stirrups (pair)	ANS001
Leg troughs (pair)	LGT001
Phlebotomy arm (single / pair)	PAS001 PAP001
Adjustable arm supports (pair)	AAS001
Options to order	
Additional +3" width	A3WB
Breather hole with vinyl wedge to plug hole	BHW001
Additional top layer of anatomic memory foam	AMF001
Alternative electronic footswitch	AEF001
Base variant to allow hoist access	BVH001

Three Section Couches

For a wider range of uses.

Multi-discipline three-section models suit a variety of medical, clinical and minor surgical applications. Electronic variable height is included as standard, with selected couches available with a hydraulic lift alternative.



Retractable wheel base with foot lever control allows couch to be moved/positioned with ease.

Attachments to order:



Features:

- Gas assisted adjustable backrest and footrest sections
- Retractable wheelbase with foot lever control
- Anti-microbial seamless 2" vinyl upholstery
- Steel frame with white anti-microbial protective coating.

Dimensions:

Length: 1900mm Height adjustment: 480-890mm Safe working load: 225kg

Vinyl Colours:

- 1 Burnt Red
- 2 Rose
- 3 Mulberry
- 4 Lilac
- 5 Navy
- 6 Royal
- 7 Lido
- 8 Aqua
- 9 Ocean
- 10 Citrus
- 11 Sage
- 12 Grey
- 13 Dove
- 14 Ginger
- 15 Buttercup
- 16 Taupe
- 17 Barley
- 18 Stone
- 19 White
- 20 Black

Options to order:



Three Section Couches

Ordering:

Add the chosen colour code to the end of the product code. For example: **CC3MD/25H6** for Royal Three Section Couch - Hydraulic

Product	Code
Three Section Couch - Electronic	CC3MD/25E
Three Section Couch - Hydraulic	CC3MD/25H
Attachments to order	
Hygiene roll holder	HRH001
Safety side rails (pair)	SSR001
Safety side pads (pair)	SSP001
Head cushion support	HCS001
IV pole	IVP001
Ankle stirrups (pair)	ANS001
Leg troughs (pair)	LGT001
Phlebotomy arm (single / pair)	PAS001 PAP001
Adjustable arm supports (pair)	AAS001
Options to order	
Additional +3" width	A3WB
Breather hole with vinyl wedge to plug hole	BHW001
Additional top layer of anatomic memory foam	AMF001
Alternative electronic footswitch	AEF001
Base variant to allow hoist access	BVH001

Paediatric Couches

For a wider range of uses.

Paediatric treatment couches are configured to smaller dimensions in electric or hydraulic models with a variety of attachment options for the medical treatment of children. The two section couches benefit from an adjustable back rest and a choice of height control operations.



Retractable wheel base with foot lever control allows couch to be moved/positioned with ease.

Attachments to order:



*CC2PA/24R only

Features:

- Retractable wheelbase with foot lever control
- Anti-microbial seamless 2" vinyl upholstery
- Steel frame with white anti-microbial protective coating.

Dimensions:

Length: 1400mm Height adjustment: 480-890mm Safe working load: 195kg

Vinyl Colours:

- 1 Burnt Red
- 2 Rose
- 3 Mulberry
- 4 Lilac
- 5 Navy
- 6 Royal
- 7 Lido
- 8 Aqua
- 9 Ocean
- 10 Citrus
- 11 Sage
- 12 Grey
- 13 Dove
- 14 Ginger
- 15 Buttercup
- 16 Taupe
- 17 Barley
- 18 Stone
- 19 White
- 20 Black

Options to order:



Paediatric Couches

Ordering:

Add the chosen colour code to the end of the product code. For example: **CC2PA/24RH14** for Ginger Paediatric Couch - Hydraulic

Product	Code
Paediatric Couch - Electronic	CC2PA/24RE
Paediatric Couch - Hydraulic	CC2PA/24RH
Attachments to order	
Hygiene roll holder	HRH001
Safety side rails (pair)	SSR001
Safety side pads (pair)	SSP001
Head cushion support	HCS001
IV pole	IVP001
Ankle stirrups (pair)	ANS001
Leg troughs (pair)	LGT001
Phlebotomy arm (single / pair)	PAS001 PAP001
Adjustable arm supports (pair)	AAS001
Options to order	
Breather hole with vinyl wedge to plug hole	BHW001
Additional top layer of anatomic memory foam	AMF001
Alternative electronic footswitch	AEF001
Base variant to allow hoist access	BVH001

Bariatric Couches

Multi-discipline two-section and three-section wider models with increased safe working load to accommodate bariatric patients in a variety of medical, clinical and minor surgical applications. Electronic variable height comes as standard, with an optional hydraulic lift alternative available on selected couches.



Attachments to order:



*CC2MD/28 & CC3MD/28 only

Features:

- Gas assisted adjustable backrest section
- Retractable wheelbase with foot lever control
- Anti-microbial seamless 2" vinyl upholstery
- Steel frame with white anti-microbial protective coating.

Dimensions:

Length: 1900mm Height adjustment: 480-890mm Safe working load: 250kg

Vinyl Colours:

- 1 Burnt Red
- 2 Rose
- 3 Mulberry
- 4 Lilac
- 5 Navy
- 6 Royal
- 7 Lido
- 8 Aqua
- 9 Ocean
- 10 Citrus
- 11 Sage
- 12 Grey
- 13 Dove
- 14 Ginger
- 15 Buttercup
- 16 Taupe
- 17 Barley
- 18 Stone
- 19 White
- 20 Black

Options to order:



Bariatric Couches

Ordering:

Add the chosen colour code to the end of the product code. For example: **CC2MD/28H5** for Navy Two Section Bariatric Couch - 28" Wide - Hydraulic

Product	Code
Two Section Bariatric Couch - 28" Wide - Electronic	CC2MD/28E
Two Section Bariatric Couch - 28" Wide - Hydraulic	CC2MD/28H
Two Section Bariatric Couch - 32" Wide - Electronic	CC2MD/32E
Two Section Bariatric Couch - 32" Wide - Hydraulic	CC2MD/32H
Three Section Bariatric Couch - Electronic	CC3MD/28E
Three Section Bariatric Couch - Hydraulic	CC3MD/28H
Attachments to order	
Hygiene roll holder	HRH001
Safety side rails (pair)	SSR001
Safety side pads (pair)	SSP001
Head cushion support	HCS001
IV pole	IVP001
Ankle stirrups (pair) (CC2MD/28 & CC3MD/28 only)	ANS001
Leg troughs (pair) (CC2MD/28 & CC3MD/28 only)	LGT001
Phlebotomy arm (single / pair)	PAS001 PAP001
Adjustable arm supports (pair)	AAS001
Options to order	
Breather hole with vinyl wedge to plug hole	BHW001
Additional top layer of anatomic memory foam	AMF001
Alternative electronic footswitch	AEF001
Base variant to allow hoist access	BVH001

Obstetrics/Gynaecology Couches

Our specialist couches have been designed for gynaecology and obstetric procedures with a maximum tilt of 25° on the end seat section and gas-assisted, hydraulically adjusted section controls. Our popular range includes both electronic and static height models.

Options:



Twin steel pillar for added rigidity

Scissor frame gives vertical elevation and stability

Attachments to order:



Obstetrics/Gynaecology Couches

Features:

- Retractable wheel base with foot lever control
- Anti-microbial seamless 2" vinyl upholstery
- Steel frame with white anti-microbial protective coating
- **Electronic Model:**
Electronic height operation by base foot rail
- **Hydraulic Model:**
Hydraulic height operation by pneumatic foot pump

Dimensions:

Length: 2030mm Height adjustment: 480 - 860mm Safe working load: 225kg

Vinyl Colours:

- | | | | | |
|-------------|-----------|------------|-----------|--------------|
| 1 Burnt Red | 2 Rose | 3 Mulberry | 4 Lilac | 5 Navy |
| 6 Royal | 7 Lido | 8 Aqua | 9 Ocean | 10 Citrus |
| 11 Sage | 12 Grey | 13 Dove | 14 Ginger | 15 Buttercup |
| 16 Taupe | 17 Barley | 18 Stone | 19 White | 20 Black |

Options to order:



AMF001

Ordering:

Add the chosen colour code to the end of the product code. For example: **CC2GN/25E4** for Lilac Gynaecology Couch - Electronic

Product	Code
Gynaecology Couch - Electronic	CC2GN/25E
Gynaecology Couch - Hydraulic	CC2GN/25H
Gynaecology Scissor Couch - Electronic	CC2SC/25E
Gynaecology Scissor Couch - Hydraulic	CC2SC/25H
Attachments to order	
Adjustable leg supports, metal tray	SUPPLIED
Hygiene roll holder	HRH001
Options to order	
Breather hole with vinyl wedge to plug hole	BHW001
Additional top layer of anatomic memory foam	AMF001

Obstetrics/Gynaecology Chairs

Our gynaecology chairs have been specially designed for all gynaecology, colposcopy and obstetric procedures, using insights from healthcare professionals.

Options:



Options to order:



Obstetrics/Gynaecology Chairs

Ordering:

Add the chosen colour code to the end of the product code. For example: **CH3GN/2413** for Dove Static Gynaecology Chair

Product	Code
Static Gynaecology Chair	CH3GN/24
Gynaecology Chair - Electronic	CH3GN/24E
Gynaecology Chair - Hydraulic	CH3GN/24H
Attachments to order	
Adjustable leg supports, arm rests, head rest, metal tray	SUPPLIED
Options to order	
Additional top layer of anatomic memory foam	AMF001

Features:

- Anti-microbial 2" vinyl upholstery
- Steel frame with white anti-microbial protective coating
- Manually adjusted leg section can be easily removed.

Static Chair:

- Static fixed height at 34"
- Gas assisted adjustable back section.

Adjustable Chair:

- Variable height
- Gas assisted adjustable back and mid sections
- Seat section has a tilt of 20°
- Retractable wheel base with foot lever control.
- **Electronic Model:**
Electronic height operation by foot console
- **Hydraulic Model:**
Hydraulic height operation by pneumatic foot pump

Dimensions:

Static Chair: Height: 860mm **Safe working load:** 195kg

Adjustable Chair: Height adjustment: 580 - 860mm **Safe working load:** 225kg

Vinyl Colours:

- | | | | | |
|--------------------|------------------|-------------------|------------------|---------------------|
| 1 Burnt Red | 2 Rose | 3 Mulberry | 4 Lilac | 5 Navy |
| 6 Royal | 7 Lido | 8 Aqua | 9 Ocean | 10 Citrus |
| 11 Sage | 12 Grey | 13 Dove | 14 Ginger | 15 Buttercup |
| 16 Taupe | 17 Barley | 18 Stone | 19 White | 20 Black |

Multi-Purpose Clinic Chairs

Versatile and practical clinical treatment chairs for surgeries, hospitals or private practice, with anti-microbial steel frame and vinyl upholstery for infection control and a choice of electronic or hydraulic operation for easy adjustment.

Options:



Retractable wheel base allows chair to be moved and positioned with ease

Attachments to order:



Options to order:



Features:

- Seat section has a fixed tilt of 10° (CH3MP/24, CH3MP/28 only)
- Gas assisted adjustable back and foot sections
- Retractable wheelbase with foot lever control.

Dimensions:

Standard: Length: 1780mm Height adjustment: 480-840mm Safe working load: 225kg
Wide: Length: 1780mm Height adjustment: 480-840mm Safe working load: 225kg
Extra Wide: Length: 1820mm Height adjustment: 480-890mm
 Safe working load: 250kg

Vinyl Colours:

- 1 Burnt Red
- 2 Rose
- 3 Mulberry
- 4 Lilac
- 5 Navy
- 6 Royal
- 7 Lido
- 8 Aqua
- 9 Ocean
- 10 Citrus
- 11 Sage
- 12 Grey
- 13 Dove
- 14 Ginger
- 15 Buttercup
- 16 Taupe
- 17 Barley
- 18 Stone
- 19 White
- 20 Black

Multi-Purpose Clinic Chairs

Ordering:

Add the chosen colour code to the end of the product code.
 For example: **CH3MP/24E17** for Barley Clinic Chair - Electronic

Product	Code
Clinic Chair - Electronic	CH3MP/24E
Clinic Chair - Hydraulic	CH3MP/24H
Wide Clinic Chair - Electronic	CH3MP/28E
Wide Clinic Chair - Hydraulic	CH3MP/28H
Bariatric Treatment Chair - Electronic	CH3MP/30
Attachments to order	
Safety side rails (pair)	SSR001
Safety side pads (pair)	SSP001
Head cushion support	HCS001
IV pole	IVP001
Phlebotomy arm (single / pair)	PAS001 PAP001
Adjustable arm supports (pair)	AAS001
Options to order	
Additional top layer of anatomic memory foam	AMF001
Alternative electric footswitch	AEF001
Base variant to allow hoist access	BVH001

Paediatric Multi-Purpose Clinic Chairs

Designed for professionals. Adapted for children.

Our compact paediatric chairs are configured to smaller dimensions for the medical treatment of children, with a variety of attachments to allow for clinical and surgical procedures, minor ops and phlebotomy. The antimicrobial steel frame and vinyl upholstery promotes hygiene while the choice of handheld remote control or pneumatic foot pump offers versatility.

Options:



Retractable wheel base allows chair to be moved and positioned with ease

Attachments to order:



*JOS4160 only

Features:

- Compact Chairs configured to smaller proportions for children
- Seat section has a fixed tilt of 10°
- Gas assisted adjustable back and foot sections
- Retractable wheelbase with foot lever control.

Dimensions:

Length: 1450mm Height adjustment: 480-840mm Safe working load: 225kg

Vinyl Colours:

- 1 Burnt Red
- 2 Rose
- 3 Mulberry
- 4 Lilac
- 5 Navy
- 6 Royal
- 7 Lido
- 8 Aqua
- 9 Ocean
- 10 Citrus
- 11 Sage
- 12 Grey
- 13 Dove
- 14 Ginger
- 15 Buttercup
- 16 Taupe
- 17 Barley
- 18 Stone
- 19 White
- 20 Black

Options to order:



Paediatric Multi-Purpose Clinic Chairs

Ordering:

Add the chosen colour code to the end of the product code.
For example: **CH3PA/21E19** for White Paediatric Childrens Clinic Chair - Electronic

Product	Code
Paediatric Childrens Clinic Chair - Electronic	CH3PA/21E
Paediatric Childrens Clinic Chair - Hydraulic	CH3PA/21H
Paediatric Childrens Phlebotomy Chair - Electronic	CH3PB/21E
Paediatric Childrens Phlebotomy Chair - Hydraulic	CH3PB/21H
Attachments to order	
Safety side rails (pair) (CH3PA/21 only)	SSR001
Safety side pads (pair) (CH3PA/21 only)	SSP001
Head cushion support	HCS001
IV pole	IVP001
Phlebotomy arm (single / pair) (CH3PA/21 only)	PAS001 PAP001
Adjustable arm supports (pair)	AAS001
Options to order	
Additional top layer of anatomic memory foam	AMF001
Alternative electric footswitch	AEF001
Base variant to allow hoist access	BVH001

Phlebotomy Chairs

Our blood sampling phlebotomy clinic chairs, with a choice of single or double arm rests, are suitable for hospitals, surgeries or private practices. The seat section has a fixed tilt of 10° (CH3PB/24, CH3PB/28 only) with a choice of electronic or hydraulic height operation controls.

Options:



Retractable wheel base allows chair to be moved and positioned with ease

Attachments to order:



*CH3PB/28, CH3PB/30 only

Options to order:



Features:

- Gas assisted adjustable back and foot sections
- Retractable wheelbase with foot lever control
- Anti-microbial seamless 2" vinyl upholstery
- Steel frame with white anti-microbial protective coating.

Dimensions:

Standard: Length: 1780mm Height adjustment: 480-840mm Safe working load: 225kg
Wide: Length: 1780mm Height adjustment: 480-840mm Safe working load: 225kg
Extra Wide: Length: 1820mm Height adjustment: 480-890mm Safe working load: 250kg

Vinyl Colours:

- 1 Burnt Red
- 2 Rose
- 3 Mulberry
- 4 Lilac
- 5 Navy
- 6 Royal
- 7 Lido
- 8 Aqua
- 9 Ocean
- 10 Citrus
- 11 Sage
- 12 Grey
- 13 Dove
- 14 Ginger
- 15 Buttercup
- 16 Taupe
- 17 Barley
- 18 Stone
- 19 White
- 20 Black

Phlebotomy Chairs

Ordering:

Add the chosen colour code to the end of the product code. For example: **CH3PB/24E11** for Sage Phlebotomy Chair - Electronic

Product	Code
Phlebotomy Chair - Electronic	CH3PB/24E
Phlebotomy Chair - Hydraulic	CH3PB/24H
Wide Phlebotomy Chair - Electronic	CH3PB/28E
Wide Phlebotomy Chair - Hydraulic	CH3PB/28H
Bariatric Phlebotomy Chair - Electronic	CH3PB/30
Attachments to order	
Phlebotomy arm rests (pair) (CH3PB/24 only)	SUPPLIED
Phlebotomy arm rest (single) (CH3PB/28, CH3PB/30 only)	SUPPLIED
Head cushion support	HCS001
IV pole	IVP001
Additional phlebotomy arm rest (CH3PB/28, CH3PB/30 only)	PAS001
Adjustable arm supports (pair)	AAS001
Options to order	
Additional top layer of anatomic memory foam	AMF001
Alternative electric footswitch	AEF001
Base variant to allow hoist access	BVH001

Vena Vario Phlebotomy Chairs

Combine comfort and seamless functionality.

Easy operation makes it possible to act quickly in emergency situations, with the patient able to be moved safely and quickly into the Trendelenburg position. The Vario Levant has 2 x 'Vena® fixed' with PU arm pads, integrated hydraulic foot pump operation and features a synchronously adjustable back and leg rest allowing the operator to place the patient in the optimum position when carrying out procedures.

Options:



Vario Phlebotomy Chair

Vario Levant Phlebotomy Chair

Features:

- Available in Stamskin® and PUxx Material options
- Wipe clean hygienic design.
- Robust steel base, white (RAL 9016) with 4 feet
- Protective glides
- Medical device in accordance with MDR 2017/745 (class I).

Vario Levant:

- Hydraulic height adjustment 540-740mm
- Hand controlled synchronic back and leg rest operation.

Stamskin® Colours:

- | | | | |
|---|---|---|---|
|  K85 Anthracite |  K05 Black |  K294 Blush |  K170 Burgundy Red |
|  K189 Capri Blue |  K89 Dark Blue |  K196 Dutch Orange |  K292 Forest Green |
|  K295 Mouse Grey |  K190 Ochre Yellow |  K394 Pink Orchid |  K258 Taupe |

Ordering:

Add the chosen colour code to the end of the product code. For example: **PC2001/MAK29** for Dark Blue Vena Vario Phlebotomy Chair.

Product	Code
Vena Vario Phlebotomy Chair	PC2001/MA
Vena Vario Levant Phlebotomy Chair	PCH2001



Ultrasound Scanning/Medical Procedure Chairs

Precise patient positioning.

Choose from electronic height operation with foot console or hydraulic height operation with pneumatic foot pump to achieve optimal patient positioning. The long backrest is ideal for ultrasound examinations, while the compact base allows close access to the patient. The steel frame and vinyl upholstery both benefit from an anti-microbial protective coating for maximum hygiene.



Provides precise patient positioning with raised seat tilt facility.

Attachments to order:



Options to order:



Ultrasound Scanning/Medical Procedure Chairs

Features:

- Variable height
- Gas assisted adjustable back and mid sections
- Manually adjusted leg section
- Seat section has hydraulic tilt of 20°
- Retractable wheelbase with foot lever control.

Dimensions:

Length: 1860mm Height adjustment: 580-860mm Safe working load: 225kg

Vinyl Colours:

- | | | | | |
|-------------|-----------|------------|-----------|--------------|
| 1 Burnt Red | 2 Rose | 3 Mulberry | 4 Lilac | 5 Navy |
| 6 Royal | 7 Lido | 8 Aqua | 9 Ocean | 10 Citrus |
| 11 Sage | 12 Grey | 13 Dove | 14 Ginger | 15 Buttercup |
| 16 Taupe | 17 Barley | 18 Stone | 19 White | 20 Black |

Ordering:

Add the chosen colour code to the end of the product code. For example: **CH3PR/24E16** for Taupe Ultrasound / Scanning / Clinic Chair - Electronic

Product	Code
Ultrasound / Scanning / Clinic Chair - Electronic	CH3PR/24E
Ultrasound / Scanning / Clinic Chair - Hydraulic	CH3PR/24H
Attachments to order	
Head cushion support	HCS001
Leg support troughs (pair)	LGT001
Adjustable arm supports (pair)	AAS001
Options to order	
Additional top layer of anatomic memory foam	AMF001

Mammography Chair

Accessible solutions for all patients.

Our Mammography Chair is designed for patients who need to sit down during examinations. This may include women prone to fainting, wheelchair users or patients with diseases affecting the musculoskeletal system. Its compound upholstery can be washed with spirits at concentrations of 70%–85% alcohol and withstands chlorine disinfection containing active chlorine at 1000–1200 ppm.



Features:

- Smooth-running wheels, push bar and foot brake for quick and secure patient positioning for the examination device
- Optimal patient positioning with electric height adjustment
- Adjustable height (510-810mm) and backrest
- Foldable armrests for transferring patients.

Dimensions:

Seat: Width (mm) 470, Depth (mm) 470

Backrest: Width (mm) 400, Height (mm) 400

Armrest Pad: Width (mm) 90, Length (mm) 260

Frame Size: 550 x 550mm

Dolce Colours:



Ordering:

Add the chosen colour code to the end of the product code.

For example: **700324** for Poppy Mammography Chair.

Product	Code
Mammography Chair	70032



“We just all love the chair here at the department. It is so easy to manoeuvre and really helps us overcome the huge patient flow of many different patients.”

We occasionally also use the chair for shorter transfers e.g., bringing the patient from the waiting room to the examination.

We are looking into the possibility of getting a chair in each of our other examination rooms because the need to handle elderly and larger patients is getting more and more of an issue for the staff.

The flexibility of the Mammography chair means we are able to examine both larger and smaller patients whilst ensuring the correct working posture.

The chair contributes to handling heavy, elderly or disabled patients faster and safer than without.

We feel the patients who use the chair for examinations are more comfortable and relaxed in what can be an uncomfortable and stressful situation.”

Vickie Sanderson

*Specialist Radiographer & Lead Mammographer
Queen Elizabeth II Hospital*