

Gynaecology Chair - Static

Product code: CH3GN/24

Agile Medical



Vinyl colours:




Attachments note:




NEW DESIGN Padded arm rests attach to sides of chair. Can be adjusted back / down when not required. Arm rests; leg supports; head rest & metal tray - all included with chair.

Specification:



Footprint 25" x 33" (63cm x 84cm)



Features:

- Chair for Gynaecology, Colposcopy, Obstetric procedures
- Static fixed height
- Gas assisted adjustable back section
- Manually adjusted leg section can be easily removed
- Anti-microbial 2" vinyl upholstery
- Steel frame with white anti-microbial protective coating



Attachments to order:

SUPPLIED AS STANDARD

Adjustable leg supports, arm rests, head rest, metal tray

Options to order:

AMF001

Additional top layer of anatomic memory foam



Manufacture and delivery:

- **FREE Delivery** mainland UK by pallet courier
- **Personal installation** service to order
- **Production lead time** approx 20 working days

For advice and ordering:

call +44 (0)1472 266 190

email: sales@agilemedical.co.uk

Gynaecology Chair - Adjustable

Agile Medical

Product code: **CH3GN/24E** (Electronic) / **CH3GN/24H** (Hydraulic)



Leg section
easily removed

Vinyl colours:




Attachments note:




NEW DESIGN Padded arm rests attach to sides of chair. Can be adjusted back / down when not required. Arm rests; leg supports; head rest & metal tray - all included with chair.

Specification:



Footprint 25" x 33" (63cm x 84cm)



Features:

- Chair for Gynaecology, Colposcopy, Obstetric procedures
- Static fixed height
- Gas assisted adjustable back section
- Manually adjusted leg section can be easily removed
- Anti-microbial 2" vinyl upholstery
- Steel frame with white anti-microbial protective coating

Electronic Model:

- Electronic height operation by foot console

Hydraulic Model:

- Hydraulic height operation by pneumatic foot pump.



UK Made. CE marked



Anti-microbial vinyl & frame

Attachments to order:

SUPPLIED AS STANDARD

Adjustable leg supports, arm rests, head rest, metal tray

Options to order:

AMF001

Additional top layer of anatomic memory foam



AMF001

Manufacture and delivery:

- **FREE Delivery** mainland UK by pallet courier
- **Personal installation** service to order
- **Production lead time** approx 20 working days

For advice and ordering:

call +44 (0)1472 266 190

email: sales@agilemedical.co.uk