

Obstetrics / Gynaecology Couch

Product code: CC2GN/25

Agile Medical



Vinyl colours:




Attachments note:




Gynae couch features perennial cut out.

Multi-adjustable leg supports and metal tray
- **included with couch**

Specification:



Footprint 25" x 65" (63cm x 165cm)



Features:

- Gynae & Obstetric Couch
- Static height
- Tubular frame gives optimum strength and stability
- Manually adjusted section controls
- Sections adjust to give a completely flat couch
- End-seat section has a max. tilt of 25°
- Anti-microbial seamless 2" vinyl upholstery
- Steel frame with white anti-microbial protective coating
- Supplied flat packed.



UK Made. CE marked



Anti-microbial vinyl & frame

Attachments to order:

HRH001

Hygiene roll holder



HRH001

Options to order:

BHW001

Breather hole with vinyl wedge to plug hole

AMF001

Additional top layer of anatomic memory foam



BHW001



AMF001

Manufacture and delivery:

- **FREE Delivery** mainland UK by pallet courier
- **Personal installation** service to order
- **Production lead time** approx 15-20 working days

For advice and ordering:

call +44 (0)1472 266 190

email: sales@agilemedical.co.uk