

Eco-Space Infusion Stand

Material: alu / stainless steel

ECO12112

PRODUCT DESCRIPTION

Field of application	Normal Care
Material	Stainless steel
Material foot base	Aluminum, ultramarine blue
Colour of plastic parts	Blue
Bottle Holder	Stainless steel, for screw-in
Number of hooks	4
Arrangement of the hooks	90° to each other
Type of height adjustment	One-hand height adjustment
Maximum length/height	2,220
Minimum length/height	1,360
Inner tube diameter	Ø 18 mm
Outer tube diameter	Ø 25 mm
Length of outer tube	1,000
Castors	Mono castors
Castor diameter	Ø 75 mm
Storage	Space-saving stackable
Space-saving	61%
Drip vessel sterilizable	Yes
Total load capacity (kg)	9.5
Load capacity per system (kg)	6
Load capacity per hook (kg)	2



PACKAGING UNIT: UNIT

Width Box1	0,59 m
Height Box1	0,25 m
Length Box1	0,64 m
Weight Box1	2,8 kg
Width Box2	0,2 m
Height Box2	0,1 m
Length Box2	1,18 m
Weight Box2	1,75 kg



Noise-reduced
high quality, low-noise castors



DIN ISO 15375
Bottle holders according to standards

ECO12112 Eco-Space Infusion Stand



Safety

Low risk of pinching when retracting and extending the IV-pole thanks to the built-in safety distance

No risk of injury thanks to the rounded shape of the bottle hooks

The geometry of the hooks prevents the infusion bottle from falling out

