

Anaesthetic Trolleys

“Our Anaesthetic trolley provides the option to anaesthetise in theatre, removing any unnecessary patient movement and improving safety. The standardised layout and content of our trolleys make it easy for practitioners to transfer between theatres and hospitals and we have eliminated a large amount of waste and out of date products by removing the majority of stock from our Anaesthetic rooms.”

Adam Dobson

*Lead Anaesthetist,
Manchester University Hospitals NHSFT*



Manchester University
NHS Foundation Trust



*Accessory to be ordered separately

Anaesthetic Trolleys



Efficient equipment wherever you need it.

Used daily by hundreds of anaesthetic practitioners throughout the NHS, our award-winning anaesthetic trolleys reduce staff downtime and prevent 'dead stock' build up by mobilising the storage of anaesthetic equipment. This compact, space-optimising solution can be configured to your specific requirements.

Multiple options of tilt bin available, as indicated over the page. Please specify your preference along with your trolley of choice.

Options:

Narrow Trolley



LX Trolley



Double Trolley



Anaesthetic Trolleys

Features:

- 10-year warranty
- Epoxy coated metal construction increases durability.
- One piece, fluid retaining ABS top for easy cleaning.
- Protective corner profiles and base bumper
- Integrated push handle
- Slide pull out extension worksurface.
- Multiple locking options for securing internal contents.

Colours:



Dimensions:

Narrow Trolley:

Width (mm) 520, Depth (mm) 540, Height (mm) 946

LX Trolley:

Width (mm) 690, Depth (mm) 540, Height (mm) 1010

Double Trolley:

Width (mm) 1125, Depth (mm) 545, Height (mm) 1010



Tilt Bin Options:

Half-size
(for Narrow
Trolley)

Small

Large

Extra large

Ordering:

Product	Drawers	Col	Code
Narrow Anaesthetic Trolley			NC34ANE1
			NC34ANE2
LX Anaesthetic Trolley			LX34ANE1
			LX34ANE2
			LX34ANE3
Double Anaesthetic Trolley			DC34ANE1
			DC34ANE2
Accessories			
Magnetic label holders			SUPPLIED
Drawer divider set (2 supplied, additional available)			SUPPLIED
Anaesthesia label dispenser			AMDSP-HC10
Tilt bin organiser (half size for Narrow Trolley)			TBS312
Tilt bin organiser (small)			TBS612
Tilt bin organiser (large)			TBL611
Tilt bin organiser (extra large)			TBX612
Sharps bin holder			SBH100
Catheter holder			CHS114
1 box glove holder			GDS119
Gel holder			GMH100
Security clamp			DRS200
Waste container with lid			WCS213

For our full range of trolley accessories see pages **108-111**.