

Bariatric Couches

Multi-discipline two-section and three-section wider models with increased safe working load to accommodate bariatric patients in a variety of medical, clinical and minor surgical applications. Electronic variable height comes as standard, with an optional hydraulic lift alternative available on selected couches.



Two Section - 28"



Two Section - 32"



Three Section - 28"

Attachments to order:



HRH001



SSR001



SSP001



HCS001



IVP001



ANS001*



LGT001*



PAS001/PAP001



AAS001

*CC2MD/28 & CC3MD/28 only

Bariatric Couches

Features:

- Gas assisted adjustable backrest section
- Retractable wheelbase with foot lever control
- Anti-microbial seamless 2" vinyl upholstery
- Steel frame with white anti-microbial protective coating.

Dimensions:

Length: 1900mm **Height adjustment:** 480-890mm **Safe working load:** 250kg

Vinyl Colours:

| | | | | |
|--------------------|------------------|-------------------|------------------|---------------------|
| 1 Burnt Red | 2 Rose | 3 Mulberry | 4 Lilac | 5 Navy |
| 6 Royal | 7 Lido | 8 Aqua | 9 Ocean | 10 Citrus |
| 11 Sage | 12 Grey | 13 Dove | 14 Ginger | 15 Buttercup |
| 16 Taupe | 17 Barley | 18 Stone | 19 White | 20 Black |

Options to order:



Ordering:

Add the chosen colour code to the end of the product code.

For example: **CC2MD/28H5** for Navy Two Section Bariatric Couch - 28" Wide - Hydraulic

| Product | Code |
|---|-------------------------------|
| Two Section Bariatric Couch - 28" Wide - Electronic | CC2MD/28E |
| Two Section Bariatric Couch - 28" Wide - Hydraulic | CC2MD/28H |
| Two Section Bariatric Couch - 32" Wide - Electronic | CC2MD/32E |
| Two Section Bariatric Couch - 32" Wide - Hydraulic | CC2MD/32H |
| Three Section Bariatric Couch - Electronic | CC3MD/28E |
| Three Section Bariatric Couch - Hydraulic | CC3MD/28H |
| Attachments to order | |
| Hygiene roll holder | HRH001 |
| Safety side rails (pair) | SSR001 |
| Safety side pads (pair) | SSP001 |
| Head cushion support | HCS001 |
| IV pole | IVP001 |
| Ankle stirrups (pair) (CC2MD/28 & CC3MD/28 only) | ANS001 |
| Leg troughs (pair) (CC2MD/28 & CC3MD/28 only) | LGT001 |
| Phlebotomy arm (single / pair) | PAS001 PAP001 |
| Adjustable arm supports (pair) | AAS001 |
| Options to order | |
| Breather hole with vinyl wedge to plug hole | BHW001 |
| Additional top layer of anatomic memory foam | AMF001 |
| Alternative electronic footswitch | AEF001 |
| Base variant to allow hoist access | BVH001 |