

Bowl Stands

Manufactured from grade 304 stainless steel, our bowl stands and kickabouts exceed infection control requirements for within sterile environments. All parts can easily be cleaned and sanitised.

Options:

Single Bowl Stand



350mm

Double Bowl Stand



700mm

Features:

- Stainless steel frame (304 grade)
- 1 or 2 removable, stainless steel basin
- 5x 75mm castors (2x braked) – stainless steel base.

Dimensions:

Single Bowl Stand:

Width (mm) 350, Depth (mm) 350, Height (mm) 800

Double Bowl Stand:

Width (mm) 700, Depth (mm) 350, Height (mm) 800

Ordering:

Product	Code
Single Bowl Stand	AM/BSD-79
Double Bowl Stand	AM/BSD-80