



Increased efficiency with reduced footprint.

The most space-efficient, general-purpose trolley on the market. Designed and manufactured in the UK, Carebuddy trolleys from Agile Medical provide the optimum solution for mobilising storage of medical supplies throughout healthcare facilities, using a tray support system with positive stops.



Features:

- Awarded Clinell® material compatibility
- Robust welded steel frame with durable epoxy coating
- One piece, fluid retaining ABS top for easy cleaning
- Protective corner bumpers.

Dimensions:


Width (mm) 664, Depth (mm) 458, Height (mm) 950

Storage Configuration:

5x Medium Polycarbonate Tray
 Width (mm) 600, Depth (mm) 400, Height (mm) 110


**Various sizes, colours and configurations possible
 (contact us today to configure your Carebuddy)**

clinell® approved



**30%
 more
 efficient
 storage***

*Compared to market alternatives



Ordering:

| Product | Col | Code |
|------------------------------|-----|--------------------|
| Carebuddy Trolley - White | | CBOB10WT050 |
| Carebuddy Trolley - Blue | | CBOB10BL050 |
| Carebuddy Trolley - Green | | CBOB10GN050 |
| Carebuddy Trolley - Orange | | CBOB10OR050 |
| Carebuddy Trolley - Red | | CBOB10RD050 |
| Carebuddy Trolley - Yellow | | CBOB10YL050 |
| Trolley Accessories | | |
| Gel holder | | GMH100 |
| Catheter holder | | CHS114 |
| Oxygen tank holder (CD size) | | OTS114 |
| Glove dispenser | | GDST19 |
| Sharps bin holder | | SBH100 |
| Waste container with lid | | WCS213 |
| IV pole | | IVP100 |

For our full range of trolley accessories see pages **108-111**.

Tray accessories:

