

Half Section Trolleys - Standard

For all your essentials.

Our Half Section Trolleys provide storage for essential equipment in a lightweight, compact and robust mobile trolley, which fits neatly under standard worktops.



Features:

- Welded, steel frame with nylon protective coating
- Basket/tray runners with positive stops
- One piece ABS top
- 4x swivel TPR castors with protective buffers.

Dimensions:

Width (mm) 350, Depth (mm) 445, Height (mm) 840

Ordering:

Product	Code
Half Section Trolley - 4 x Wire Baskets	304100A
Half Section Trolley - 1 x Poly Tray, 3 x ABS Trays	304100B
Half Section Trolley - 3 x ABS Baskets, 1 x Wire Basket	304100C
Accessories	
Sharps bin holder	306155
Foldable side frame	307045
Vertical catheter holder	301106
Horizontal catheter holder	301107
Notes holder	301310
400mm push handle	601014
Apron roll dispenser	306171

Half Section Trolleys - Dublin

Features:

- Welded, steel frame with nylon protective coating
- Basket/tray runners with positive stops
- One piece ABS top
- 4x swivel TPR castors with protective buffers.

Dimensions:

Width (mm) 350, Depth (mm) 445, Height (mm) 840

Ordering:

Product	Code
Half Section Dublin Trolley - 5 x Poly Baskets	304114A
Half Section Dublin Trolley - 5 x ABS Trays	304114B
Half Section Dublin Trolley - 5 x ABS Baskets	304114C
Accessories	
Sharps bin holder	306155
Foldable side frame	307045
Vertical catheter holder	301106
Horizontal catheter holder	301107
Notes holder	301310
400mm push handle	601014
Apron roll dispenser	306171

