

Heavy Duty Infusion Stand

For the most challenging environments

Our most resilient solution for intensive care medicine with a soft drop system that dampens and controls the lowering of the inner tube. Each of the 4 stainless steel hook can carry loads of up to 5kg with 5 x 75mm twin castors providing a secure and stable base.



Features:

- Stainless steel tubes according to DIN EN 10217-7
- Material ultrasonically and eddy current tested (according to EN 10246-9/3)
- Height adjustable 1430 – 2260mm
- Total load capacity 35kg.

Ordering:

Product	Code
Heavy Duty Infusion Stand	PC1172954
Accessories	
Pump Holder	Z2N06292
Shelf	K2S7043A
Push handle	Z2N06982



Pump holder



Shelf



Push handle