

Height Adjustable Chart Trolley

Support where you need it.

A lightweight, compact trolley for maintaining patient charts with 360mm of hydraulic height adjustment to easily accommodate all users. Choose from our range of accessories to customise your workstation to your specific needs.



Features:

- Aluminium trolley with acrylic work surface
- Secure chart clips and integral pen tray
- 25kg load capacity
- 4 x 100mm lockable castors (2 x braked).

Dimensions:

Trolley Height:

775 - 1135mm

Worksurface:

Width (mm) 620, Height (mm) 560

Chart Board:

Width (mm) 450, Height (mm) 450

Wheelbase Diameter:

746mm

Ordering:

Product	Code
Height Adjustable Chart Trolley	L6-200-CH