

Infection Control Trolley

Features:

- Easy clean surface
- Lockable swivel castors and push handle
- Constructed using high pressure laminate panels.
- 2 x shelves
- 3 box glove holder
- Hand disinfectant dispenser
- Drip bin

Dimensions:

Width (mm) 420, Depth (mm) 540, Height (mm) 850
Height including back panel: 1435mm

Ordering:

Product	Code
Infection Control Trolley	PEC05111C1

A mobile solution to stop the spread.

Our closed storage cupboard is a simple yet effective solution to infection control. With standard accessories including a glove box holder and hand gel dispenser.

