

Instrument Trolley with Drawers

Features:

- Durable, impact-resistant steel construction with powder coating
- 1x 76mm, 1x 152mm drawers with full extension runners
- 4x 75mm castors (2x braked).

Dimensions:

Width (mm) 500, Depth (mm) 500, Height (mm) 900

Ordering:

Product	Code
Instrument Trolley with Drawers	IT3502SW
Accessories	
ABS Top	IT5050T

Maximising even the smallest of spaces.

The Slimline Instrument Trolley is almost half the width of our standard instrument trolleys and is equipped with 2 storage drawers and storage shelves with side rails, making it ideal for hospitals and healthcare centres. The optional, moulded ABS top has fluid retaining edges for easy cleaning.

