

Instrument Trolleys - One Shelf

Low maintenance, high strength.

Our Tall Instrument Trolley is seamlessly welded from 304 grade stainless steel for strength. Rounded corners replace sharp edges to minimize risk of injury, while lockable swivel castors provide stability.

Options:

900mm High Trolley

1000mm High Trolley



Features:

- Easy clean stainless steel
- Fixed shelf with down turned sides.
- 4x 75mm swivel non marking castors (2x braked).

Dimensions:

900mm High:

Width (mm) 900, Depth (mm) 450, Height (mm) 900

1000mm High:

Width (mm) 900, Depth (mm) 450, Height (mm) 1000

Ordering:

Product	Code
Instrument Trolley - One Shelf - 900mm High	SIT0994
Instrument Trolley - One Shelf - 1000mm High	SIT1094