

Instrument Trolleys

Seamless strength with a softer edge.

Our Instrument Trolleys are constructed using seamlessly welded strong stainless steel for ease of cleaning with rounded corners and thermoplastic buffers for safety.

Options:

Instrument Trolley (450mm width)



Instrument Trolley (610mm width)



Instrument Trolley (750mm width)



Instrument Trolleys

Features:

- Fixed shelves with down turned sides.
- Maximum capacity per shelf 40kg
- 4x 75mm swivel non-marking, anti-static castors.

Dimensions:

450mm Width Trolley:

Width (mm) 450, Depth (mm) 450, Height (mm) 840

610mm Width Trolley:

Width (mm) 610, Depth (mm) 450, Height (mm) 840

750mm Width Trolley:

Width (mm) 750, Depth (mm) 450, Height (mm) 840

910mm Width Trolley:

Width (mm) 910, Depth (mm) 450, Height (mm) 840

Ordering:

Product	Code
Instrument Trolley - 450x450mm	SS4545F/B
Instrument Trolley - 610x450mm	SS6145F/B
Instrument Trolley - 750x450mm	SS7545F/B
Instrument Trolley - 910x450mm	SS9145F/B

Instrument Trolley (910mm width)

