

Large Procedure Trolley

An innovative new addition.

Features:

- Durable epoxy-coated metal construction with 10-year warranty
- Easy clean, one piece, fluid retaining ABS top.
- Protective corner profiles and base bumper
- Slide pull out extension work surface.

Dimensions:

Width (mm) 930, Depth (mm) 525, Height (mm) 1020

Ordering:

Product	Col	Code
Large Procedure Trolley		VI34PRO1
Accessories		
Magnetic label holders		SUPPLIED
IV pole		IVP100
Sharps bin holder		SBH100
Catheter holder		CHS114
1 box glove holder		GDS119
Gel holder		GMH100
Security clamp		DRS200
Waste container with lid		WCS213
Multi-storage basket		MSL100
File holder		MSL100

For our full range of trolley accessories see pages **108-111**.

The newest product in our award-winning procedure trolley range, the Large Procedure Trolley provides the solution to ensuring secure storage of long procedural equipment, most commonly laparoscopes. It measures 930mm wide with integrated push handle and full width drawers that have multiple locking options for secure storage.

