

Light Duty Coloured Infusion Stand

Features:

- 4 plastic hooks (load 2kg per hook)
- Screw height adjustment 1350 – 2150mm
- Total load capacity 9.5kg
- 5x 50mm castors.

Colours:



Ordering:

Product	Col	Code
Light Duty Coloured Infusion Stand		I-N51121
		I-N61121
		I-N71121
		I-N81121
		I-N91121
Accessories		
Shelf		K2S7043A
Push handle		Z2N06982



Shelf Push handle

A brighter option for paediatric care.

Our child friendly Coloured Infusion stands may look fun, but they provide serious support for Paediatric professionals. The patented sealing cap prevents ingress of hygienically sensitive liquids, while the robust 3-layer nickel and chrome construction provides durability with high corrosion protection and excellent salt spray test resistance.

