

Light Duty Infusion Stand

Safety and security as standard.

Manufactured to the highest standards, using ultrasonically and eddy current tested material (according to EN 10246-9/3) and DIN EN 10217-7 certified stainless steel tubes. 5 x 50mm castors provide stability, while the rounded shape of the bottle hooks prevents risk of injury.



Features:

- Stainless steel tubes according to DIN EN 10217-7
- Material ultrasonically and eddy current tested (according to EN 10246-9/3)
- Height adjustable 1430 – 2260mm
- Total load capacity 9.5kg.

Ordering:

Product	Code
Light Duty Infusion Stand	NS17229
Accessories	
Shelf	K2S7043A
Push handle	Z2N06982



Shelf

Push handle