

Mayo Tables

Features:

- Stainless steel Frame (304 Grade)
- Stainless steel top with perimetrical rim
- Height adjustable with choice of hydraulic / foot pedal or manual adjustment
- 400 x 600mm or 500 x 700mm top
- 4x 75mm castors (2x braked).

Height adjustment:

Manual Height Adjustment:

800 - 1200mm

Hydraulic Height Adjustment:

900 - 1200mm

Ordering:

Product	Code
Mayo Table - Manual Height Adjustment 400 x 600mm top	HMTD45
Mayo Table - Manual Height Adjustment 500 x 700mm top	HMTD46
Mayo Table - Hydraulic Height Adjustment 400 x 600mm top	HMTD49
Mayo Table - Hydraulic Height Adjustment 500 x 700mm top	HMTD50

Easy to use, mobile Mayo Tables providing a place for sterile instruments and supplies used during surgery. Fully stainless-steel welded construction providing compatibility with all standard hospital cleaning products.



Manual Height
Adjustment

Hydraulic Height
Adjustment