

Medical Grade Keyboards



Pathoban Medical Keyboard

Features:

- IP68 water and dust protection
- 100% silicone body (latex free)
- USB wired version with waterproof silicone USB cap
- Real scissor type keys with excellent tactile responsiveness
- Plug and play – no drivers needed
- Cleaning lock key and cleaning reminder function
- 0.5mm raised surfaces make it easy to clean.

Ordering:

Product	Code
Pathoban Medical Keyboard	SK10801
SterileFLAT Antibacterial Standard Keyboard	SF0902V4
SterileFLAT Antibacterial Backlit Keyboard	SF20-BL
K10 Medical Grade Keyboard	MCW-K10



SterileFLAT Antibacterial Standard Keyboard

Features:

- USB Keyboard. Layout: UK
- Material: Nano Silver Silicon membrane
- Compatible with: Win98/Me/2000/NT/XP/Vista/Windows 7, 8, 8.1, 10, 11
- Membrane: Latex Free
- Numeric Keypad: Yes
- Mounting: Desktop and Medical Cart
- Approval Certification: CE, FCC, Approved RoHS compliant
- Cable Length 1.5m
- Multimedia function: NO
- Soft touch tactile key
- Warranty: 1 Year.



SterileFLAT Antibacterial Backlit Keyboard

Features:

- USB Keyboard. Layout: UK
- 2 light level blue led backlighting
- Material: Nano Silver Silicon membrane
- Compatible with: Win98/Me/2000/NT/XP/Vista/Windows 7/8/8.1/10/11
- Membrane: Latex Free
- Numeric Keypad: Yes
- Mounting: Desktop and Medical Cart
- Approval Certification: CE, FCC, Approved RoHS compliant
- Cable Length 1.5m
- Multimedia function: NO
- Soft touch tactile keys
- Warranty: 1 Year.



K10 Medical Grade Keyboard

Features:

- USB Keyboard. Layout: UK
- Complete protection against dust and liquids IP68
- Washable, sterilised in just seconds, no seams or crevices
- Resistant against all common disinfectants and detergents
- Low force and precise feedback key technology
- Rugged shock protected chassis
- Low noise typing
- Compact size, ergonomically low profile design
- Damage protected key design
- Integrated high resolution capacitive touchpad
- CleanLock-function to avoid unintended inputs (on option)
- Warranty: 1 Year.
- Weight: 1kg
- No. of Keys: 105.

Medical Grade Mice



SterileMOUSE Antibacterial Laser Mouse

Features:

- Dimensions: 105 x 61 x 38 mm
- Weight: 100g.



Pathoban Sealed Washable Mouse

Features:

- Fully sealed
- USB
- IP68 water and dust protection, chemical resistant
- Durable silicone rubber material
- Capacitive touch sensor in middle area, 3 standard mouse left key and right key, middle key
- Colour: White.



Medigenic Scroll Mouse

Features:

- 5 buttons (left, middle, right, scroll up and scroll down)
- Dimensions: 113 x 67 x 37mm
- Weight: 106g
- Cable length: 2m
- Colour: Light blue.



IP68 Washable Antibacterial Scroll Wheel Mouse

Features:

- Nano silver antibacterial technology
- IP68 complete protection against dust and liquids
- Washable with soapy water or in dishwasher
- Can be used with medical gloves
- Resistant against most common detergents and disinfection materials
- Sterilised in seconds, no need for expensive sterilisation methods
- 1000 DPI
- Plug and play USB
- Low noise.

Ordering:

Product	Code
SterileMOUSE Antibacterial Laser Mouse	SF08-14
Pathoban Sealed Washable Mouse	ASMTI
Medigenic Scroll Mouse	MSM-03
IP68 Washable Antibacterial Scroll Wheel Mouse	STWM042