

Medical Stools

Specialist seating solutions.

Specially developed for use in medical environments such as hospital examination rooms, emergency departments and operating rooms. The optional Ergo Shape seat with recess edge helps prevent the stagnation of blood circulation to the thighs.

Options:

Round Seat



Ergo Seat



Optional foot operated height adjustment











Features:

- Height Adjustment 550 – 740mm
- Seat ø 400 mm
- Aluminium base ø 540 mm
- 65mm brake unloaded soft castors.

Seat Options:

Round or Ergo Seat

Stamskin® Colours:

-  K85 Anthracite
-  K05 Black
-  K294 Blush
-  K170 Burgundy Red
-  K189 Capri Blue
-  K89 Dark Blue
-  K196 Dutch Orange
-  K292 Forest Green
-  K295 Mouse Grey
-  K190 Ochre Yellow
-  K394 Pink Orchid
-  K258 Taupe

Ordering:

Add the chosen colour code to the end of the product code. For example: **ME6301K196** for Dutch Orange Medical Stool - Ergo Seat.

| Product | Code |
|---------------------------------|---------------|
| Medical Stool - Round Seat | MR6101 |
| Medical Stool - Ergo Seat | ME6301 |
| Accessories | |
| Foot Operated Height Adjustment | FA6000 |