

Medication Trolleys



Customisable, compact storage solutions.

Our award-winning medication trolleys, with integrated push handles, are an innovative space-saving solution for transporting equipment during ward rounds. Multiple locking options ensure security of drugs and internal contents with the opportunity to choose between 1, 2 or 3 drawers with medication bin compartments and extendable work surfaces.

Options:

Narrow Trolley



LX Trolley



Medication Trolleys

Features:

- 10-year warranty
- Durable epoxy-coated metal construction
- Easy clean, one piece, fluid retaining ABS top.
- Protective corner profiles and base bumper.

Colours:



Dimensions:

Narrow Trolley:

Width (mm) 520, Depth (mm) 540, Height (mm) 946

LX Trolley:

Width (mm) 760, Depth (mm) 490, Height (mm) 1010

Ordering:

Product	Drawers	Col	Code
Narrow Medication Trolley	1 x Medication Drawer		NCM3411
			NCM3412
	2 x Medication Drawers		NCM3421
			NCM3422
	3 x Medication Drawers		NCM3431
			NCM3432
LX Medication Trolley	1 x Medication Drawer		LXM3411
			LXM3412
			LXM3413
	2 x Medication Drawers		LXM3421
			LXM3422
			LXM3423
	3 x Medication Drawers		LXM3431
			LXM3432
			LXM3433
Accessories			
Medbin drawers and medbins with dividers			SUPPLIED
Security clamp			DRS200
IV pole			IVP100
Sharps bin holder			SBH100
Glove dispenser			GDST19
Gel holder			GMH100

For our full range of trolley accessories see pages **108-111**.