

Multi-Purpose Clinic Chairs

Versatile and practical clinical treatment chairs for surgeries, hospitals or private practice, with anti-microbial steel frame and vinyl upholstery for infection control and a choice of electronic or hydraulic operation for easy adjustment.

Options:



Standard 24"

Wide 28"

Extra Wide 30"

Retractable wheel base allows chair to be moved and positioned with ease

Attachments to order:



Options to order:



Multi-Purpose Clinic Chairs

Features:

- Seat section has a fixed tilt of 10° (CH3MP/24, CH3MP/28 only)
- Gas assisted adjustable back and foot sections
- Retractable wheelbase with foot lever control.

Dimensions:

Standard: Length: 1780mm Height adjustment: 480-840mm Safe working load: 225kg

Wide: Length: 1780mm Height adjustment: 480-840mm Safe working load: 225kg

Extra Wide: Length: 1820mm Height adjustment: 480-890mm

Safe working load: 250kg

Vinyl Colours:



Ordering:

Add the chosen colour code to the end of the product code.

For example: **CH3MP/24E17** for Barley Clinic Chair - Electronic

Product	Code
Clinic Chair - Electronic	CH3MP/24E
Clinic Chair - Hydraulic	CH3MP/24H
Wide Clinic Chair - Electronic	CH3MP/28E
Wide Clinic Chair - Hydraulic	CH3MP/28H
Bariatric Treatment Chair - Electronic	CH3MP/30
Attachments to order	
Safety side rails (pair)	SSR001
Safety side pads (pair)	SSP001
Head cushion support	HCS001
IV pole	IVP001
Phlebotomy arm (single / pair)	PAS001 PAP001
Adjustable arm supports (pair)	AAS001
Options to order	
Additional top layer of anatomic memory foam	AMF001
Alternative electric footswitch	AEF001
Base variant to allow hoist access	BVH001