

# Obstetrics/Gynaecology Chairs

Our gynaecology chairs have been specially designed for all gynaecology, colposcopy and obstetric procedures, using insights from healthcare professionals.

## Options:



Options to order:



# Obstetrics/Gynaecology Chairs

## Features:

- Anti-microbial 2" vinyl upholstery
- Steel frame with white anti-microbial protective coating
- Manually adjusted leg section can be easily removed.

## Static Chair:

- Static fixed height at 34"
- Gas assisted adjustable back section.

## Adjustable Chair:

- Variable height
- Gas assisted adjustable back and mid sections
- Seat section has a tilt of 20°
- Retractable wheel base with foot lever control.
- **Electronic Model:**  
Electronic height operation by foot console
- **Hydraulic Model:**  
Hydraulic height operation by pneumatic foot pump

## Dimensions:

**Static Chair: Height:** 860mm **Safe working load:** 195kg

**Adjustable Chair: Height adjustment:** 580 - 860mm **Safe working load:** 225kg

## Vinyl Colours:



## Ordering:

Add the chosen colour code to the end of the product code.

For example: **CH3GN/2413** for Dove Static Gynaecology Chair

| Product   | Code             |
|---|------------------|
| Static Gynaecology Chair                                  | <b>CH3GN/24</b>  |
| Gynaecology Chair - Electronic                            | <b>CH3GN/24E</b> |
| Gynaecology Chair - Hydraulic                             | <b>CH3GN/24H</b> |
| Attachments to order                                      |                  |
| Adjustable leg supports, arm rests, head rest, metal tray | <b>SUPPLIED</b>  |
| Options to order  |                  |
| Additional top layer of anatomic memory foam              | <b>AMF001</b>    |