

Obstetrics/Gynaecology Couches

Our specialist couches have been designed for gynaecology and obstetric procedures with a maximum tilt of 25° on the end seat section and gas-assisted, hydraulically adjusted section controls. Our popular range includes both electronic and static height models.

Options:



CC2GN/25

Twin steel pillar for added rigidity



CC2SC/25

Scissor frame gives vertical elevation and stability

Attachments to order:



HRH001

Obstetrics/Gynaecology Couches

Features:

- Retractable wheel base with foot lever control
- Anti-microbial seamless 2" vinyl upholstery
- Steel frame with white anti-microbial protective coating
- **Electronic Model:**
Electronic height operation by base foot rail
- **Hydraulic Model:**
Hydraulic height operation by pneumatic foot pump

Dimensions:

Length: 2030mm **Height adjustment:** 480 - 860mm **Safe working load:** 225kg

Vinyl Colours:



Options to order:



AMF001

Ordering:

Add the chosen colour code to the end of the product code.

For example: **CC2GN/25E4** for Lilac Gynaecology Couch - Electronic

Product	Code
Gynaecology Couch - Electronic	CC2GN/25E
Gynaecology Couch - Hydraulic	CC2GN/25H
Gynaecology Scissor Couch - Electronic	CC2SC/25E
Gynaecology Scissor Couch - Hydraulic	CC2SC/25H
Attachments to order	
Adjustable leg supports, metal tray	SUPPLIED
Hygiene roll holder	HRH001
Options to order	
Breather hole with vinyl wedge to plug hole	BHW001
Additional top layer of anatomic memory foam	AMF001