

# One Section Trolleys - Standard

## Broader options for essential storage.

Choose from poly trays or wire baskets to store essential equipment in our lightweight, compact and robust mobile trolley which fits discreetly under standard worktops.

**306100A**



**306100B**



**306100C**



## Features:

- Welded, steel frame with nylon protective coating
- Basket/tray runners with positive stops
- One piece ABS top
- 4x swivel TPR castors with protective buffers.

## Dimensions:

Width (mm) 650, Depth (mm) 445, Height (mm) 840

## Ordering:

Product	Code
One Section Trolley - 4 x ABS Baskets	<b>306100A</b>
One Section Trolley - 3 x Poly Trays, 1 x Wire Basket	<b>306100B</b>
One Section Trolley - 3 x Poly Trays, 1 x ABS Basket	<b>306100C</b>
Accessories	
Sharps bin holder	<b>306155</b>
Foldable side frame	<b>307045</b>
Vertical catheter holder	<b>301106</b>
Horizontal catheter holder	<b>301107</b>
Notes holder	<b>301310</b>
400mm push handle	<b>601014</b>
Apron roll dispenser	<b>306171</b>

# One Section Trolleys - Dublin

## Features:

- Welded, steel frame with nylon protective coating
- Basket/tray runners with positive stops
- One piece ABS top
- 4x swivel TPR castors with protective buffers.

## Dimensions:

Width (mm) 650, Depth (mm) 445, Height (mm) 840

## Ordering:

Product	Code
One Section Dublin Trolley - 5 x Poly Trays	<b>306114A</b>
One Section Dublin Trolley - 5 x ABS Trays	<b>306114B</b>
One Section Dublin Trolley - 5 x ABS Baskets	<b>306114C</b>
Accessories	
Sharps bin holder	<b>306155</b>
Foldable side frame	<b>307045</b>
Vertical catheter holder	<b>301106</b>
Horizontal catheter holder	<b>301107</b>
Notes holder	<b>301310</b>
400mm push handle	<b>601014</b>
Apron roll dispenser	<b>306171</b>

