

Paediatric Multi-Purpose Clinic Chairs

Designed for professionals. Adapted for children.

Our compact paediatric chairs are configured to smaller dimensions for the medical treatment of children, with a variety of attachments to allow for clinical and surgical procedures, minor ops and phlebotomy. The antimicrobial steel frame and vinyl upholstery promotes hygiene while the choice of handheld remote control or pneumatic foot pump offers versatility.

Options:



Clinic Chair



Phlebotomy Chair

Retractable wheel base allows chair to be moved and positioned with ease

Attachments to order:



*JOS4160 only

Paediatric Multi-Purpose Clinic Chairs

Features:

- Compact Chairs configured to smaller proportions for children
- Seat section has a fixed tilt of 10°
- Gas assisted adjustable back and foot sections
- Retractable wheelbase with foot lever control.

Dimensions:

Length: 1450mm **Height adjustment:** 480-840mm **Safe working load:** 225kg

Vinyl Colours:



Options to order:



Ordering:

Add the chosen colour code to the end of the product code.

For example: **CH3PA/21E19** for White Paediatric Childrens Clinic Chair - Electronic

| Product | Code |
|--|-------------------------------|
| Paediatric Childrens Clinic Chair - Electronic | CH3PA/21E |
| Paediatric Childrens Clinic Chair - Hydraulic | CH3PA/21H |
| Paediatric Childrens Phlebotomy Chair - Electronic | CH3PB/21E |
| Paediatric Childrens Phlebotomy Chair - Hydraulic | CH3PB/21H |
| Attachments to order | |
| Safety side rails (pair) (CH3PA/21 only) | SSR001 |
| Safety side pads (pair) (CH3PA/21 only) | SSP001 |
| Head cushion support | HCS001 |
| IV pole | IVP001 |
| Phlebotomy arm (single / pair) (CH3PA/21 only) | PAS001 PAP001 |
| Adjustable arm supports (pair) | AAS001 |
| Options to order | |
| Additional top layer of anatomic memory foam | AMF001 |
| Alternative electric footswitch | AEF001 |
| Base variant to allow hoist access | BVH001 |