

Paediatric Resuscitation Trolleys



Reliable, durable and... lovable!

Our paediatric-themed trolleys with cute characters add a child-friendly feel to our reliable Resuscitation Trolley. Durable epoxy-coated metal construction, secure locking options and easy clean ABS tops make these trolleys as practical as they are pretty.

Options:

Narrow Trolley

LX Trolley

Double Trolley



520mm



690mm



1125mm

Paediatric Resuscitation Trolleys

Features:

- 10-year warranty
- Protective corner profiles and base bumper
- Integrated push handle
- Slide pull out extension worksurface.

Themes:



Garden Theme



Sea Theme

Dimensions:

Narrow Trolley:

Width (mm) 520, Depth (mm) 540, Height (mm) 946

LX Trolley:

Width (mm) 690, Depth (mm) 540, Height (mm) 1010

Double Trolley:

Width (mm) 1125, Depth (mm) 545, Height (mm) 1010

Ordering:

| Product | Code |
|---|-------------------|
| Narrow Paediatric Resuscitation Trolley - Garden | NC34GAR/R2 |
| Narrow Paediatric Resuscitation Trolley - Sea | NC34SEA/R2 |
| LX Paediatric Resuscitation Trolley - Garden | LX34GAR/R2 |
| LX Paediatric Resuscitation Trolley - Sea | LX34SEA/R2 |
| Double Paediatric Resuscitation Trolley - Garden | DC34GAR/R2 |
| Double Paediatric Resuscitation Trolley - Sea | DC34SEA/R2 |
| Accessories | |
| Breakaway seals | SUPPLIED |
| Magnetic label holders | SUPPLIED |
| Drawer divider set (2 supplied, additional available) | SUPPLIED |
| Oxygen tank holder (CD size) | OTS114 |
| Stainless steel defibrillator shelf | DSS433 |
| IV pole | IVP100 |
| Sharps bin holder | SBH100 |
| Catheter holder | CHST14 |
| T box glove holder | GDST19 |
| Gel holder | GMH100 |
| Security clamp | DRS200 |
| Waste container with lid | WCS213 |

For our full range of trolley accessories see pages **108-111**.