

# Paediatric Trolleys



**Serious storage with a fun side.**

Our paediatric trolleys, with integrated push handle, have all of the benefits of the Agile Medical trolleys but are finished with attractive pictures designed to help children relax in hospital. Multiple locking options provide secure storage for internal contents and the one-piece fluid-retaining ABS top is easily cleanable.

**Options:**

**Narrow Trolley**



520mm

**LX Trolley**



690mm

# Paediatric Trolleys

## Features:

- Durable epoxy-coated metal construction with 10-year warranty
- Protective corner profiles and base bumper
- Slide pull out extension worksurface.

## Themes:



## Colours:



## Dimensions:

### Narrow Trolley:

Width (mm) 520, Depth (mm) 540, Height (mm) 946

### LX Trolley:

Width (mm) 760, Depth (mm) 490, Height (mm) 1010

## Ordering:

Product	Col	Code
Narrow Paediatric Trolley - Garden		<b>NC34GAR1</b>
		<b>NC34GAR2</b>
Narrow Paediatric Trolley - Sea		<b>NC34SEA1</b>
		<b>NC34SEA2</b>
LX Paediatric Trolley - Garden		<b>LX34GAR1</b>
		<b>LX34GAR2</b>
LX Paediatric Trolley - Sea		<b>LX34SEA1</b>
		<b>LX34SEA2</b>
Accessories		
Magnetic label holders		<b>SUPPLIED</b>
Drawer divider set (2 supplied, additional available)		<b>SUPPLIED</b>
IV pole		<b>IVP100</b>
Sharps bin holder		<b>SBH100</b>
Catheter holder		<b>CHST14</b>
1 box glove holder		<b>GDST19</b>
Gel holder		<b>GMH100</b>
Security clamp		<b>DRS200</b>
Waste container with lid		<b>WCS213</b>
Multi-storage basket		<b>MSL100</b>
File holder		<b>MSL100</b>

For our full range of trolley accessories see pages **108-111**.