

# Phlebotomy Chairs

Our blood sampling phlebotomy clinic chairs, with a choice of single or double arm rests, are suitable for hospitals, surgeries or private practices. The seat section has a fixed tilt of 10° (CH3PB/24, CH3PB/28 only) with a choice of electronic or hydraulic height operation controls.

## Options:



Standard 24"

Wide 28"

Extra Wide 30"

Retractable wheel base allows chair to be moved and positioned with ease

## Attachments to order:



\*CH3PB/28, CH3PB/30 only

## Options to order:



# Phlebotomy Chairs

## Features:

- Gas assisted adjustable back and foot sections
- Retractable wheelbase with foot lever control
- Anti-microbial seamless 2" vinyl upholstery
- Steel frame with white anti-microbial protective coating.

## Dimensions:

**Standard:** Length: 1780mm Height adjustment: 480-840mm Safe working load: 225kg

**Wide:** Length: 1780mm Height adjustment: 480-840mm Safe working load: 225kg

**Extra Wide:** Length: 1820mm Height adjustment: 480-890mm

**Safe working load:** 250kg

## Vinyl Colours:



## Ordering:

Add the chosen colour code to the end of the product code.

For example: **CH3PB/24E11** for Sage Phlebotomy Chair - Electronic

Product	Code
Phlebotomy Chair - Electronic	<b>CH3PB/24E</b>
Phlebotomy Chair - Hydraulic	<b>CH3PB/24H</b>
Wide Phlebotomy Chair - Electronic	<b>CH3PB/28E</b>
Wide Phlebotomy Chair - Hydraulic	<b>CH3PB/28H</b>
Bariatric Phlebotomy Chair - Electronic	<b>CH3PB/30</b>
Attachments to order	
Phlebotomy arm rests (pair) (CH3PB/24 only)	<b>SUPPLIED</b>
Phlebotomy arm rest (single) (CH3PB/28, CH3PB/30 only)	<b>SUPPLIED</b>
Head cushion support	<b>HCS001</b>
IV pole	<b>IVP001</b>
Additional phlebotomy arm rest (CH3PB/28, CH3PB/30 only)	<b>PAS001</b>
Adjustable arm supports (pair)	<b>AAS001</b>
Options to order	
Additional top layer of anatomic memory foam	<b>AMF001</b>
Alternative electric footswitch	<b>AEF001</b>
Base variant to allow hoist access	<b>BVH001</b>