

# Phlebotomy Trolley

## Ready to roll IT solutions.

Our Phlebotomy Cart is a solution to the ever-increasing need to bring IT to the point of care in hospitals. Built-in multi storage options accommodate nearly everyone's needs.



## Features:

- Central key locking system
- 3-year warranty
- Easy clean ABS top with anti-bacterial coating
- TPR medical castors, 2 braked.

## Dimensions:

Width (mm) 520, Depth (mm) 540, Height (mm) 946

## Ordering:

Product	Code
Phlebotomy Trolley	<b>PNC-01</b>
Accessories	
ABS bins and dividers	<b>SUPPLIED</b>
Sharps bin holder	<b>SBH100</b>
1 box glove holder	<b>GDS119</b>
Gel holder	<b>GMH100</b>
Waste container with lid	<b>WCS213</b>
Multi-storage basket	<b>MSL100</b>
File holder	<b>FHS531</b>
Bulk storage bins (large)	<b>BSB711</b>

For our full range of trolley accessories see pages **108-111**.

**Our needs-based approach enables us to create a ready to roll solution, which minimises downtime and reduces demand on internal resources.**

Leeds Teaching Hospitals NHS Trust benefitted from 53 Agile Medical Phlebotomy Trolleys, complete with laptop, printer and professionally mounted hand sanitisers.



\*Accessory to be ordered separately