

# Procedure Trolleys



A winner for all departments.

Suitable for applications hospital-wide, the Procedure Trolley combines the benefit of compact mobile storage with excellent manoeuvrability and integrated push handle. Moving your equipment storage to a trolley with lockable drawers improves security and hygiene by limiting unwanted access and potential build-up of dust and contaminants.

## Options:



# Procedure Trolleys

## Features:

- Durable epoxy-coated metal construction with 10-year warranty
- Easy clean, one piece, fluid retaining ABS top.
- Protective corner profiles and base bumper
- Slide pull out extension worksurface.

## Colours:



## Dimensions:

### Narrow Trolley:

Width (mm) 520, Depth (mm) 540, Height (mm) 946

### LX Trolley:

4 Drawers: Width (mm) 690, Depth (mm) 540, Height (mm) 930

5 Drawers: Width (mm) 690, Depth (mm) 540, Height (mm) 1010

### Double Trolley:

Width (mm) 1125, Depth (mm) 545, Height (mm) 1010

## Ordering:

Product	Drawers	Col	Code
Narrow Procedure Trolley			<b>NC34051</b>
			<b>NC34052</b>
LX Procedure Trolley	4 (1x3", 2x6", 1x9")		<b>LX31PRO1</b>
			<b>LX31PRO2</b>
	5 (2x3", 2x6", 1x9")		<b>LX34PRO1</b>
			<b>LX34PRO2</b>
Double Procedure Trolley			<b>DC34091</b>
			<b>DC34092</b>
Accessories			
Magnetic label holders			<b>SUPPLIED</b>
Drawer divider set (2 supplied, additional available)			<b>SUPPLIED</b>
IV pole			<b>IVP100</b>
Sharps bin holder			<b>SBH100</b>
Catheter holder			<b>CHST14</b>
1 box glove holder			<b>GDST19</b>
Gel holder			<b>GMH100</b>
Security clamp			<b>DRS200</b>
Waste container with lid			<b>WCS213</b>
Multi-storage basket			<b>MSL100</b>
File holder			<b>MSL100</b>

For our full range of trolley accessories see pages **108-111**.