

Resuscitation Trolleys

Leading NHS Foundation Trust, Guy's and St Thomas looked to Agile Medical to supply 200 resuscitation trolleys, to support outstanding patient care in emergency scenarios.

Guy's and St Thomas NHS Foundation Trust staff now have immediate access to appropriate resuscitation equipment and drugs to facilitate rapid resuscitation of patients in cardiorespiratory arrest. Standardisation of the equipment used for resuscitation (including defibrillators and emergency suction equipment) throughout the trust has helped to ensure that every nurse is comfortable with the resuscitation trolley.



Guy's and St Thomas'
NHS Foundation Trust



*Accessory to be ordered separately

LX Resuscitation Trolleys



Innovative design for ultimate reliability.

The most reliable Resuscitation Trolley in the UK, with its unique 'Zero-Risk' tamper evident locking bar, was designed alongside healthcare workers, to accommodate all the required emergency equipment in one secure trolley. Staff simply open one of the drawers to automatically break the tamper-evident seal, allowing quick and risk-free access to equipment. The durable epoxy-coated metal construction, combined with protective corner profiles, ensures full protection of both the trolley and the building.

Options:



LX Resuscitation Trolleys

Features:

- 15-year warranty
- One piece, fluid retaining ABS top for easy cleaning.
- Integrated push handle
- Slide pull out extension worksurface.
- Multiple locking options for securing internal contents.

Colours:



Dimensions:

4 Drawers:

Width (mm) 690, Depth (mm) 540, Height (mm) 930

5 Drawers:

Width (mm) 690, Depth (mm) 540, Height (mm) 1010



Ordering:

Product	Drawers	Col	Code
LX Resuscitation Trolley	4 (4x6")		LX31RES3
			LX31RES5
	5 (2x3", 2x6", 1x9")		LX34RES3
			LX34RES5
	5 (1x3", 4x6")		LX3414/R3
			LX3414/R5
Accessories			
Breakaway seals			SUPPLIED
Magnetic label holders			SUPPLIED
Drawer divider set (2 supplied, additional available)			SUPPLIED
Oxygen tank holder (CD size)			OTS114
Stainless steel defibrillator shelf			DSS433
IV pole			IVP100
Sharps bin holder			SBH100
Catheter holder			CHST14
1 box glove holder			GDST19
76mm adjustable drawer divider set			LXA076
152mm adjustable drawer divider set			LXA152
ABS trays (small)			LXT764
ABS trays (large)			LXT164

For our full range of trolley accessories see pages **108-111**.

Narrow and Double Resuscitation Trolleys



A range of options for different applications.

With the same reliability and durability as the LX Resuscitation Trolley, the Narrow and Double options provide for different applications. The Narrow Resuscitation Trolley offers a compact, space-optimising solution where space is a premium or less equipment is required to be stored, while the twin-compartment storage of the Double Resuscitation Trolley brings optimum flexibility with equipment storage and is often used to combine resus and airway equipment storage.

Options:

Narrow Trolleys

NC3414/R5



520mm

NC34RES5



520mm

Double Trolley

DC34RES5



1125mm

Narrow and Double Resuscitation Trolleys

Features:

- 15-year warranty
- One piece, fluid retaining ABS top for easy cleaning.
- Slide pull out extension worksurface.
- Multiple locking options for securing internal contents.

Colours:



Dimensions:

Narrow Trolley:

Width (mm) 520, Depth (mm) 540, Height (mm) 946

Double Trolley:

Width (mm) 1125, Depth (mm) 545, Height (mm) 1010

Ordering:

Product	Drawers	Col	Code
Narrow Resuscitation Trolley	5 (1x3", 4x6")		NC3414/R2
			NC3414/R5
	5 (2x3", 2x6", 1x9")		NC34RES2
			NC34RES5
Double Resuscitation Trolley		DC34RES2	
		DC34RES5	
Accessories			
Breakaway seals			SUPPLIED
Magnetic label holders			SUPPLIED
Drawer divider set (2 supplied, additional available)			SUPPLIED
Oxygen tank holder (CD size)			OTS114
Stainless steel defibrillator shelf			DSS433
IV pole			IVP100
Sharps bin holder			SBH100
Catheter holder			CHST14
T box glove holder			GDST19

For our full range of trolley accessories see pages **108-111**.

