

# Sackholders

**Strength, style and durability.**

Our capsule Sackholders are plastic clinical waste bins with a hands-free frame, removable body and foot pedal operation. Specifically designed and manufactured to meet the NHS demand for a plastic sackholder with absolutely no chance of rusting.

## Options:

28 Litre



42 Litre



80 Litre



# Sackholders

## Features:

- 100% durable plastic construction
- Easy clean surfaces and radius corners
- 100% corrosion resistant foot unit
- Available in all standard waste segregation colours.

## Lid/Frame Colour:



## Dimensions:

### 28 Litre:

Width (mm) 228, Depth (mm) 351, Height (mm) 671  
Weight: 6kg

### 42 Litre:

Width (mm) 444, Depth (mm) 391, Height (mm) 565  
Weight: 8kg

### 80 Litre:

Width (mm) 444, Depth (mm) 391, Height (mm) 830  
Weight: 10kg

## Ordering:

Product	Col	Code
28 Litre Sackholder - White Lid/Frame		AM28LWHWL
28 Litre Sackholder - Yellow Lid/Frame		AM28LWHYL
28 Litre Sackholder - Orange Lid/Frame		AM28LWHOL
28 Litre Sackholder - Black Lid/Frame		AM28LWHOL
28 Litre Sackholder - Tiger Stripe Lid/Frame		AM28LWHTL
42 Litre Sackholder - White Lid/Frame		AM42LWHWL
42 Litre Sackholder - Yellow Lid/Frame		AM42LWHYL
42 Litre Sackholder - Orange Lid/Frame		AM42LWHOL
42 Litre Sackholder - Black Lid/Frame		AM42LWHBL
42 Litre Sackholder - Tiger Stripe Lid/Frame		AM42LWHTL
80 Litre Sackholder - White Lid/Frame		AM80LWHWL
80 Litre Sackholder - Yellow Lid/Frame		AM80LWHYL
80 Litre Sackholder - Orange Lid/Frame		AM80LWHOL
80 Litre Sackholder - Black Lid/Frame		AM80LWHBL
80 Litre Sackholder - Tiger Stripe Lid/Frame		AM80LWHTL

All models have a white body.

## Optional Vinyl Artwork:

- No vinyl artwork
- CLINICAL WASTE ONLY
- DOMESTIC WASTE ONLY