

Sepsis Trolley

Supports delivery of SEPSIS 6.

Robust, compact trolley designed for use when treating Sepsis. Each trolley is supplied with labelled drawers to support delivery of SEPSIS 6 within 1 hour of sepsis recognition. With its unique 'Zero-Risk' tamper evident locking bar, staff simply open one of the drawers to automatically break the tamper-evident seal, allowing quick and risk-free access to equipment.



Features:

- 10-year warranty
- Epoxy coated metal construction increases durability
- One piece, fluid retaining ABS top for easy cleaning
- Protective corner profiles and base bumper
- Slide pull out extension worksurface.

Dimensions:

Width (mm) 520, Depth (mm) 540, Height (mm) 946

Ordering:

Product	Code
Sepsis Trolley	NC34SEP9
Accessories	
Oxygen tank holder (CD size)	OTS114
Stainless steel defibrillator shelf	DSS433
IV pole	IVP100
Sharps bin holder	SBH100
Catheter holder	CHST14
Glove dispenser	GDS119
Security clamp	DRS200
Waste container with lid	WCS213

For our full range of trolley accessories see pages **108-111**.