

Soiled Linen Trolleys

A quiet and durable solution for healthcare environments.

Developed in accordance with NHS infection prevention standards, the Agile Medical range of soiled linen trolleys ensure that the everyday task of collecting and transporting soiled linen is done safely, hygienically and easily thanks to 4 swivel castors and an optional foot-operated lid mechanism.

No foot operated mechanism



Foot operated mechanism



Components designed to reduce noise



ABS lid with integrated hinge increases durability

Soiled Linen Trolleys

Features:

- Lightweight, anodised aluminium body with moulded ABS base and lid
- Bag-holder with front opening for easy removal of linen bags
- Bag grip-ring with rubber corners to hold linen bags in place
- Suitable for disinfection.

Dimensions:

Single Trolley:

Width (mm) 390, Depth (mm) 410, Height (mm) 880

Double Trolley:

Width (mm) 785, Depth (mm) 410, Height (mm) 880

Triple Trolley:

Width (mm) 1180, Depth (mm) 410, Height (mm) 880

Soiled Linen Bags:

- Washable polyester soiled linen bag
- Capacity: approx. 15 kg (90 ltrs)
- Fitted with eyelets and draw cord with plastic toggle.
- Available in a range of colours to allow compliance with local, national and international guidelines.

Accessories to order:



ABS lids



Castor brake



Push handle



Soiled linen bags

Ordering:

Product	Code
Single Soiled Linen Trolley	STS348/NP
Single Soiled Linen Trolley - Foot operated mechanism	STS348
Double Soiled Linen Trolley	DTD748/NP
Double Soiled Linen Trolley - Foot operated mechanism	DTD748
Triple Soiled Linen Trolley	TTT1148/NP
Triple Soiled Linen Trolley - Foot operated mechanism	TTT1148
Accessories	
ABS lid - White	SLL340
ABS lid - Blue	SLL340B
ABS lid - Yellow	SLL340G
Castor brake	YRU75/80F
Single push handle	STS4MA
Double push handle	STS1MA
Triple push handle	STS2MA
White soiled linen bag	ASACCOBIANCOS
White soiled linen bag - Blue stripes	ASACCOBLUS
White soiled linen bag - Yellow stripes	ASACCOGIALLOS
White soiled linen bag - Red stripes	ASACCOROSSOS
White soiled linen bag - Green stripes	ASACCOVERDES