

# Stainless Steel Dressing Trolleys

Seamless strength with a softer edge.

Our Dressing Trolleys are constructed using seamlessly welded strong stainless steel for ease of cleaning with rounded corners and thermoplastic buffers for safety.

**Options:**

Dressing Trolley (450mm width)



Dressing Trolley (610mm width)



Dressing Trolley (750mm width)



# Stainless Steel Dressing Trolleys

## Features:

- Fixed shelves with down turned sides.
- Maximum capacity per shelf 40kg
- 4x 75mm swivel non-marking, anti-static castors.

## Dimensions:

### 450mm Width Trolley:

Width (mm) 450, Depth (mm) 450, Height (mm) 840

### 610mm Width Trolley:

Width (mm) 610, Depth (mm) 450, Height (mm) 840

### 750mm Width Trolley:

Width (mm) 750, Depth (mm) 450, Height (mm) 840

### 910mm Width Trolley:

Width (mm) 910, Depth (mm) 450, Height (mm) 840

## Ordering:

Product	Code
Stainless Steel Dressing Trolley - 450x450mm	<b>SS4545F/B</b>
Stainless Steel Dressing Trolley - 610x450mm	<b>SS6145F/B</b>
Stainless Steel Dressing Trolley - 750x450mm	<b>SS7545F/B</b>
Stainless Steel Dressing Trolley - 910x450mm	<b>SS9145F/B</b>

