

Sterile Case Carts - Aluminium

Hygienic distribution and collection.

Sterile case carts have been designed for the distribution of sterilised goods and the collection of dirty goods which require sterilising.



Dimensions:

Depth (mm) 710, Height (mm) 1550

Rail System Options:

Nylon Fibreglass Auto-stop Rail System – suitable for metal sterile containers placed on wire shelves and for plastic containers with lids for dirty goods.

Ordering:

Product	Width	Code
Aluminium Sterile Case Cart - 2 x Fixed Shelves	980mm	AM1622O3F
	1230mm	AM1625O10F
Aluminium Sterile Case Cart - 2 x Adjustable Shelves	980mm	AM1622O3A
	1230mm	AM1625O10A
Aluminium Sterile Case Cart - 3 x Fixed Shelves (MOQ of 5 pieces)	980mm	AM1623O3F
	1230mm	AM1626O10F
Aluminium Sterile Case Cart - 3 x Adjustable Shelves (MOQ of 5 pieces)	980mm	AM1623O3A
	1230mm	AM1626O10A
Nylon Fibreglass Auto-stop Rail System	980mm	AM3162O1
	1230mm	AM3165O4

For sterile instrument container options see page 158.