

# Sterile Case Carts

Engineered for efficiency.

Drive efficiencies in your instrument transportation process, with case carts engineered for high-capacity transport and ease of use.

Prominent visual indicators displaying cart status (optional)



Integrated push handles

Double walled, 270° wing doors keep carts open during automated washing to ensure full sterilisation

Robust construction with smooth, easy to clean surfaces. Available in steel or aluminium

Reduce transport time by maximising storage with adjustable shelves

Protective rubber bumper around the bottom edge

Strong non-twisting frame ensuring carts are suitable for towing and for transporting outside the hospital

160mm medical castors for a smooth, quiet, and stable ride

2x fixed castors to ensure directional safety when towing multiple case carts

Circumferential door seals to ensure protection of internal sterile contents

