

Sterile Case Carts - Stainless Steel

Dimensions:

Depth (mm) 710, Height (mm) 1510

Rail System Options:

Stainless Steel Rail System - Suitable for the transport of DIN-type metal containers.

Shelf Options:

Adjustable, Perforated Shelves. The shelves are perforated and embossed for optimal fluid evacuation during washing. Only compatible with stainless steel case carts.

Ordering:

Product	Width	Code
Stainless Steel Sterile Case Cart - 2 Fixed, Unperforated Shelves	1025mm	AM8622X2
	1235mm	AM8625X2
Stainless Steel Sterile Case Cart - 2 Adjustable, Perforated Shelves	1025mm	AM8622X1
	1235mm	AM8625X1
Stainless Steel Sterile Case Cart - 3 Fixed, Unperforated Shelves	1025mm	AM8623X2
	1235mm	AM8626X2
Stainless Steel Sterile Case Cart - 3 Fixed, Unperforated Shelves	1025mm	AM8623X1
	1235mm	AM8626X1
Stainless Steel Auto-stop Rail System	1025mm	AM8162X1
	1235mm	AM8165X1

For sterile instrument container options see page **158**.

Adapted for changing needs.

The increased use of alkaline products in sterilisation procedures has led to the development of a new range of case carts, entirely constructed from stainless steel.

Adjustable,
Perforated Shelves

