

Sterile Packing Table

Features:

- Manual height adjustment 741 – 1041mm
- Tested to 90kg working load.
- Can be manufactured to bespoke fabricated sizes to accommodate individual needs.

Bespoke solutions for individual needs.

Our stainless steel table provides a variable working height for the comfort of the user. The manually controlled, rise and fall tables. can be easily adjusted in height.

Ordering:

Product	Width	Code
600mm Deep Sterile Packing Table	1200mm	AMSPT1260
	1500mm	AMSPT1280
	1800mm	AMSPT1560
800mm Deep Sterile Packing Table	1200mm	AMSPT1580
	1500mm	AMSPT1860
	1800mm	AMSPT1880

